

## Large-scale full-length HTT Q19 and Q42 purification from HEK293T – 2018/04/16

### Rationale:

Previously, small scale expression of constructs TOC009-D01, D02, D04 and D06 (Q23, Q54, Q78 and Q145 respectively) was completed <https://zenodo.org/record/1239023> and showed good expression for the 2 shorter polyQ-length huntingtin constructs. Next this protocol was scaled for large scale purification of huntingtin Q19 and Q42.

### Protocol:


1250 mL HEK293T production for TOC009-F01 and 800 mL HEK293T production for F06 (HTT<sup>1-3144</sup> pBACMAM C-terminal FLAG – Q19 and Q42 respectively) – completed with baculoviral transduction protocol by Alma Seitova <https://zenodo.org/record/1194724>. Cells were harvested by centrifugation at 1500 rpm, 30 mins at 4 °C in a benchtop centrifuge and were stored at -80 °C prior to purification.

Cell pastes were thawed and resuspended in ~200 mL 50 mM Tris pH 8, 500 mM NaCl supplemented with 1 x protease inhibitors supplemented with benzonase. NB: freeze-thaw cycle for cells is sufficient for lysis. The lysate was clarified by centrifugation and then bound to 5 mL anti-FLAG resin (Sigma M2) at 4 °C with rocking for 2 hours. Resin was washed with 500 mL 50 mM Tris pH 8, 500 mM NaCl. HTT protein was eluted with ~10 mL resuspension buffer supplemented with 200 µg/mL 3xFLAG peptide. The sample was concentrated to 1 mL and run on Superose6 column. Samples were run on SDS-PAGE 4-20 % tris-glycine gel of the purification process.

Superose6 elution profile

Zoom of peaks

4-20 % Tris-Glycine SDS-PAGE



### Comments

Peak shape from elution is the typical huntingtin distribution. Yields very good and samples clean.

- F01 -> 3.5 mg/mL – 4 x 20 µL aliquots
- F06 -> 3.2 mg/mL – 4 x 20 µL aliquots