

# Purification of proteins for pulldown experiment with XIAP

## Rationale:

XIAP is known to bind to the kinase RIPK2. RIPK2 is one of the few kinase proteins that shows overlap of reactivity to some of the compounds that we've been looking at in the context of ALK2 and other type I/ type II proteins in the BMP/TGF $\beta$  families. Therefore we were interested to see if the similarity between RIPK2 and BMP/TGF $\beta$  receptors extended into them also binding to XIAP.

## Aim:

Purify the following proteins for pulldown experiments:

TGFBR1, TGFBR2, ACVR2, BMPR2, BMPR1A, BMPR1B.

## Protein Purification

### Proteins to be purified:

#### TGFBR1

MGHHHHHHSSGVDLGTENLYFQ/\*SMEDPSLDRPFISEGTTLKDLIYDMTTSGSGSGLPLLQRTIARTIVLQESIGK  
GRFGEVWRGKWRGEEVAVKIFSSREERSWFREAEIYQTVMLRHENILGFIAADNKDNGTWTQLWLVS DYHEHGS  
LFDYLNRYTVTVEGMIKLALSTASGLAHLHMEIVGTQGKPAIAHRDLKSKNILVKKNGTCCIADLGLAVRHDSATDTI  
DIAPNHRVGTKRYMAPEVLDD SINMKHFESFKRADIYAMGLVFW EIAARRCSIGGIHEDYQLPYYDLVPSDPSVEEM  
RKVVCEQKLRPNIPNRWQSCEALRVM AKIMREC WYANGAARLTALRIKKTLSQLSQQEGIKM

#### TGFBR2

MGHHHHHHSSGVDLGTENLYFQ/\*SMHNTELLPIELDTLVGKGRFAEVYKAKLKQNTSEQFETVAVKIFPYEEYAS  
WKTEKDIFSDINLKHENILQFLTAEERKTEL GKQYWLITAFHAKGNLQEYLTRHVISWEDLRKLGSS LARGIAHLHSD  
HTPCGRPKMPIVHRDLKSSNILVKNDLTCCLCD FGLSLRLDPTLSVDDLANS GGQVGTARYMAPEVLES RMNLENVE  
SFKQTDVYSMALVLWEMTSRCNAVGEVKDYEPFPGSKVREHPCVESMKDNVLRDRGRPEIPSWL NHQGIQMV  
CETLTECWDH DPEARLTAQCVAERFSELEHLDR L

#### ACVR2

MGHHHHHHSSGVDLGTENLYFQ/\*SMPLQLLEVKARGRFGCVWKAQLLNEYVAVKIFPIQDKQSWQNEYEVYSL  
PGMKHENILQFIGAEKRGTSVDVDLWLITAFHEKGSLSDFLKANVVS WNELCHIAETMARGLAYLHEDIPGLKDGH  
KPAISHRDIKSKNVLLKNNLTACIADFGLALKFEAGKSAGDTHGQVGTTRYMAPEVLEGAINFQRDAFLRIDMYAM  
GLVLWELASRCTAADGPVDEYMLPFEEEIGQHPSLEDMQEVV VHKKRVPVLRDYWQKHAGMAMLCETIEECWD  
HDAEARLSAGCVGERITQMQLRTNI

#### BMPR2

MEAAASEPSLDLDNLKLELIGRGRYGAVYKGS LDERPVAVKVFSFANRQNFINEKNIYRVPLMEHDNIARFIVGDE  
RVTADGRMEYLLVMEYYPNGSLCKYLSLHTSDWVSSCRLAHSVTRGLAYLHTELPRGDHYKPAISHRDLNSRNVLV  
KNDGTCVISDFGLSMRLTGNRLV R PGEEDNAAISEVGTIRYMAPEVLEGAVNLRDCESALKQVDMYALGLIYWEIF  
MRCTDLFPGESVPEYQMAFQTEVGNHPTFEDMQVLVSREKQRPKFPEAWKENS LAVRSLKETIEDCWDQDAEA  
RLTAQCAEERMAELMMIWERNKSVSPTV NPMSTAMQNERAHHHHHH

#### BMPR1A

MGHHHHHHSSGVDLGTENLYFQ/\*SMAFIPVGESLKDLDIQSQSSGSGSGLPLLQRTIAKQIQMV RQVVGKGRYG  
EVWVGKWRGEKVAVKVFFTEEASWFRETEIYQTVLMRHENILGFIAADIKGTGSWTQLYLITDYHENGSLYDFLK  
CATLDTRALLKLAYS AACGLCHLHTEIYGTQGKPAIAHRDLKSKNILIKKNGSCCIADLGLAVKFNSDTNEVDVPLNTR  
VGTKRYMAPEVLDES LNKNHFQPYIMADIYSFGLI WEMARRCITGGIVEEYQLPYYNMVPSDPSYEDMREVV CVK  
RLRPVSNRWNSDECLRAVLKLMSECAWHPASRLTALRIKKT LAKMVESQDVKI

## BMPR1B

MGHHHHHSSGVDLGTENLYFQ/\*SMTYIPPGESLRDLIEQSQSSGSGSGLPLLQRTIAKQIQMVKQIGKGRYGE  
VWMGKWRGEKVAVKVFFTTTEASWFRETEIYQTVLMRHENILGFIAADIKGTGSWTQLYLITDYHENGSLYDYLKS  
TTLDAKSMKLKLAYSSVSGLCHLHTEIFSTQGKPAIAHRDLKSKNILVKKNGTCCIADLGLAVKFISDTNEVDIPPNTRV  
GTKRYMPPEVLDES LNRNHFQSYIMADMYSFGLILWEVARRCVSGGIVEEYQLPYHDLVPSDPSYEDMREIVCIKKL  
RPSFPNRWSSDECLRQMGKLMTECWAHNPA SRLTALRVKKT LAKMSESQDIKL

/\* denotes Tev cleavage site


## Expression:

TGFBR1, BMPR1A, BMPR1B, TGFBR2 and ACVR2 were expressed in Sf9 insect cells by Katarzyna Kupinska at 27°C for 72 hours in glass shaker flasks.

BMPR2 was expressed in *E. coli* Rossetta cells – An initial overnight culture was grown at 37C with Kanamycin. 10ml was used to inoculate 1 L LB at 37C until OD600 1.3 (accidental overgrowth), induced with 0.4mM IPTG and subsequently grown overnight at 18C before harvesting.



## Purification:

- Thaw pellets in luke warm water
- Sonicate BMPR2A for 5 min (5 on, 10 off) on ice.
- Sonicate all other samples for 3 min (5 on, 10 off) on ice.
- Divide each sample into two centrifugation tubes and top up with 'binding buffer' (500mM NaCl, 50mM HEPES, 5mM Imidazole, 5% glycerol, pH7.5) so each tube is full (aprox 40ml per tube).
- Add 1ml PEI (5%) per tube to precipitate DNA.
- Spin at 21.5k rpm for 50 minutes.
- Incubate lysate with 2ml pre-equilibrated Ni-NTA beads at 4C for 1h. (500mM NaCl, 50mM HEPES, 5mM Imidazole, 5% glycerol, pH7.5)
- Spin down beads at 700g for 10 minutes to separate beads from lysate.
- Pour off lysate.
- Resuspend beads in 50ml binding buffer and load onto gravity flow column (collect flow through)
- Wash with 30ml wash buffer (500mM NaCl, 50mM HEPES, 30mM Imidazole, 5% glycerol, pH7.5).
- Elute in 7ml elution buffer 1 (500mM NaCl, 50mM HEPES, 50mM Imidazole, 5% glycerol, pH7.5).
- Elute in 3 x 7ml elution buffer 4 (500mM NaCl, 50mM HEPES, 250mM Imidazole, 5% glycerol, pH7.5).
- Add 1mM TCEP to all fractions.
- Run samples on an SDS PAGE gel (mix 5ul of loading dye with 15ul sample, boil for 3 minutes and load 10ul onto the gel.) and run at 160V for 50 minutes.
- Pool elution fractions E1 - E4(3)
- Add Tev protease to fractions which contain protein (Wash 2, E1, E4<sub>1</sub>)
- Incubate at 4C overnight.





*Samples corresponding to the nickle column fractions run on an SDS PAGE – FT (flowthrough combined with binding buffer wash), W (wash with wash buffer), E1 (elution with E1 elution buffer), E4<sub>(1-3)</sub> (three elution's with elution 4 buffer)*

- Concentrate down boxed samples to <5ml and run on a pre-equilibrated GF 200 column (BMPR2 and TGFBR1) or a GF 75 column (ACVR2 or TGFBR2) using standard gel filtration buffer (300mM NaCl, 50mM HEPES, 5mM Imidazole, 5% glycerol, pH7.5).
- Collect fractions and run a gel of peak.




*BMPR2 and TGFB1 SDS PAGE gel (top) shows samples from the relevant size exclusion column fractions, corresponding to the UV peak shown in the bottom panel.*




*ACVR2 and TGFB2 SDS PAGE gel (top) shows samples from the relevant size exclusion column fractions, corresponding to the UV peak shown in the bottom panel.*

- Mass spec proteins.
- Pool fractions back together and concentrate down before flash freezing in liquid nitrogen.


### BMPR2




### TGFBR1



### ACVR2



### TGFBR2



Mass spec of purified proteins.