

Institut-Hôpital neurologique de Montréal

Montreal Neurological Institute-Hospital



iPSC Protocols

Culture Reagents

Authors: Carol Xiuqing Chen, Chanshuai Han, Meghna Mathur, Nguyen-Vi Mohamed, and Cecilia Rocha

v2.0

EDDU-001-02

March 2020

Antibiotic-Antimycotic

Gibco #15240-062

 Contains the antibiotics penicillin and st amphotericin B

and st Select to go to the Table of Contents screen Select to go to the Materials section

As supplied/	100x		Solution		Select to return to this Help screen	
stock solution	1000	Columbia		Colution	date of manufacture	
Final solution (media)	iPSC culture; DA neuron, cortical neuron, and motor neuron differentiation	1x	•	Dilute 1:100	4°C	
Final solution (coating)	Laminin solution	1x	•	Dilute 1:100	4°C	





Select to navigate to the **previous screen** or to the **next screen**



Table of Contents



Introduction	Reagents
Objectives	Culture reagents: A–H
Technical and safety considerations	Culture reagents: I–P
Materials	Culture reagents: Q–Z
Culture reagents: A–H	
Culture reagents: I–P	
Culture reagents: Q–Z	









Introduction





EDDU-001-02

March 2020



Introduction

Objectives

- This document provides key information for reagents used in culture media for:
 - Peripheral blood mononuclear cell (PBMC) reprogramming
 - Induced pluripotent stem cell (iPSC) culture
 - Mouse embryonic fibroblast (MEF) culture
 - Induction of dopaminergic (DA) neural progenitor cells (NPCs), cortical NPCs, and motor NPCs
 - NPC culture
 - Differentiation of dopaminergic (DA) neurons, cortical neurons, motor neurons, cholinergic neurons, sensory neurons, GABAergic neurons, astrocytes, microglia, and oligodendrocytes
 - 3D culture of forebrain, midbrain, and cerebral organoids
- The following information is provided for each reagent as supplied and for stock and final solutions, where applicable:
 - Molecular weight (MW) of solid powders or concentration of solutions
 - Preparation instructions
 - Storage conditions
 - Important information for handling and usage









Introduction



Technical and Safety Information

- This information in this document is not intended to replace the product information or material safety data sheet (MSDS) from the manufacturer. Refer to all documentation provided by the manufacturer prior to using any reagent.
- Follow all laboratory safety requirements as outlined by McGill University Environmental Health and Safety (EHS).
- All reagents in this document are for use in cell culture and should be kept sterile. All solutions should be prepared using sterile technique.
- For reagents supplied as solid powders, reconstitute the entire vial or bottle to prepare stock solutions when appropriate.
- Final solution information is provided for the most common application(s) and is subject to change depending on the application and in the event of protocol updates. Prepare final solutions as per the most recent version of the protocol.









Materials





EDDU-001-02

March 2020



Materials



Select a culture reagent to go to its screen. Culture Reagents A-H

- Activin-A
- <u>Antibiotic-Antimycotic</u>
- Ascorbic acid (AA)
- B-27 supplement
- B-27 supplement, minus Vitamin A
- Biotin
- Bone morphogenetic protein 4 (BMP4)
- Bone morphogenetic protein 9 (BMP9)/GDF-2
- Bovine serum albumin (BSA)
- Brain-derived neurotrophic factor (BDNF)
- <u>CD200</u>
- <u>CHIR-99021</u>
- Ciliary neurotrophic factor (CNTF)
- <u>Compound E</u>
- <u>CX3CL1</u>
- <u>DAPT</u>

Institut-Hôpital neurologique de Montréal Montreal Neurological Institute-Hospital

- Dexamethasone
- Dibutyryl-cAMP (db-cAMP)

- <u>DMEM/F12</u>
- <u>DMEM</u>
- Dorsomorphin
- Dorsomorphin homologue 1 (DMH1)
- Embryonic stem-cell fetal bovine serum (FBS)
- Epidermal growth factor (EGF)
- Erythropoietin (EPO)
- Essential 8 supplement
- Essential 8 basal media
- Fetal bovine serum (FBS)
- Fibroblast growth factor-basic (FGF-b)/FGF-2
- Fibroblast growth factor 8 (FGF-8)
- <u>5-Fluoro-2'-deoxyuridine</u>
- Gentle Cell Dissociation Reagent
- Glial-derived neurotrophic factor (GDNF)
- GlutaMAX-I
- <u>Heparin</u>

🐯 McGill



Materials



Select a culture reagent to go to its screen. Culture Reagents: I-P

- Insulin
- Insulin-like growth factor 1 (IGF-1)
- Insulin-transferrin-sodium selenite supplement (ITS)
- Interleukin-3 (IL-3)
- Interleukin-6 (IL-6)
- Interleukin-34 (IL-34)
- Knockout DMEM/F12
- Knockout serum replacement
- Laminin (Invitrogen)
- Laminin (Sigma)
- L-Ascorbic acid 2-phosphate (AA2P)
- <u>LDN193189</u>
- L-Glutamine
- Lithium chloride (LiCl)
- Macrophage colony stimulating factor (M-CSF)
- <u>Matrigel growth factor reduced (GFR) basement</u> membrane matrix

- Matrigel human embryonic stem cell (hESC)-qualified matrix
- MEM nonessential amino acid (NEAA) solution
- <u>2-mercaptoethanol (β-mercaptoethanol [BME]; Gibco</u> <u>#21985023)</u>
- <u>2-mercaptoethanol (β-mercaptoethanol [BME]; Merck</u> <u>#8057400005)</u>
- <u>Mitomycin C</u>
- <u>mTeSR1 5x supplement</u>
- mTeSR1 basal media
- <u>N-2 supplement</u>
- <u>β-Nerve growth factor (β-NGF)</u>
- <u>Neurobasal (NB) media</u>
- <u>Neurotrophin-3 (NT3)</u>
- Noggin
- Penicillin-Streptomycin
- <u>Platelet-derived growth factor-AA (PDGF-AA)</u>
- Phosphate-buffered saline (PBS)
- Poly-L-ornithine (PO)
- Polyvinyl alcohol (PVA)
- Purmorphamine

Institut-Hôpital neurologique de Montréal Montreal Neurological Institute-Hospital

🐺 McGill



EDDU-001-02



Materials

Select a culture reagent to go to its screen. Culture Reagents: Q–Z

- <u>Retinoic acid (RA)</u>
- <u>SB431542</u>
- <u>Smoothened Agonist (SAG)</u>
- Sodium butyrate
- Sonic hedgehog (SHH)
- Sonic hedgehog (SHH; C24II)
- Stem cell factor (SCF)
- Stemline II hematopoietic stem cell expansion medium

McGill

- <u>StemPro Accutase Cell Dissociation Reagent</u>
- <u>Thiazovivin</u>

Institut-Hôpital neurologique de Montréal Montreal Neurological Institute-Hospital

- <u>1-Thioglycerol</u>
- Thrombopoietin (TPO)
- <u>Transforming growth factor-β3 (TGF-β3)</u>
- <u>3,3',5-Triiodo-L-thyronine (T3)</u>

- Valproic acid (VPA)
- <u>XAV939</u>
- <u>Y-27632</u>





Culture Reagents: A–H





EDDU-001-02

March 2020



Activin A



Peprotech #120-14

• Regulates cell proliferation and differentiation and promotes neuronal survival

	\sim				
	MW/Concentra	ation	Preparation instructions	Storage	Important information
As supplied	26.0 kDa		 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	12.5 µg/mL	-	 Dissolve 100 ug in 8 mL sterile ddH₂O containing 0.1% BSA Prepare 100-µL aliquots in 0.7-mL tubes 	–80°C	
Final solution (media)	Microglia differentiation	12.5 ng/mL	Dilute stock solution 1:1000	4°C	







Antibiotic-Antimycotic



Gibco #15240-062

 Contains the antibiotics penicillin and streptomycin and the antifungal agent amphotericin B

	~			
	MW/Concentration	Preparation instructions	Storage	Important information
As supplied/ stock solution	100x	Solution	–20°C for up to 1 year from date of manufacture	
Final solution (media)	iPSC culture; DA neuron, cortical neuron, and motor neuron differentiation	 Dilute stock solution 1:100 	4°C	
Final solution (coating)	Laminin solution	 Dilute stock solution 1:100 	4°C	







Ascorbic acid (AA)

Sigma #A5960

()

Institut-Hôpital neurologique de Montréal Montreal Neurological Institute-Hospital **McGill**

• An enzymatic cofactor and antioxidant

	~~				
	MW/Concentration	F	Preparation instructions	Storage	Important information
As supplied	176.12 g/mol	• F	Powder	Room temperature	May darken in storage
Stock solution	200 mM	5 • F • F	Dissolve 176.12 mg in 5 mL sterile ddH ₂ O Filter solution Prepare 200-µL aliquots in 0.7-mL tubes	–80°C	
Final	DA neuron or cortical neuron 20 differentiation	μM • [Dilute stock solution 1:1000	4°C	
solution (media)	Motor neuron differentiation, 10 midbrain organoid culture	μM • Ε	Dilute stock solution 1:2000	4°C	






B-27 supplement

Gibco #17504044

• Supplement that supports the neuronal growth and viability











B-27 supplement, minus Vitamin A



Gibco #12587010

Supplement that is ideal for the cultivation of neural progenitor and stem cells, either as
organoids in suspension or in adherent monolayer culture, without inducing
differentiation

	\sim				
	MW/Concentrat	ion	Preparation instructions	Storage	Important information
As supplied/ stock solution	50x		Solution	–20°C for 12 months from date of manufacture	Protect from light
Final solution (media)	Forebrain and midbrain organoid culture	1x	Dilute stock solution 1:50	4°C	Protect from light
	Cerebral organoid culture	0.5x	Dilute stock solution 1:100	4°C	Protect from light







Biotin



Sigma #B4639

• Essential vitamin that is important for amino acid and energy metabolism and fatty acid synthesis

	~				
	MW/Concentra	tion	Preparation instructions	Storage	Important information
As supplied	244.31 g/mc	bl	• Powder	4°C	
Stock solution	100 µg/mL		 Dissolve 10 mg in 1 mL 1N NaOH Add 100 µL of solution to 9.9 mL 1x PBS Filter solution Prepare 100-µL aliquots in 0.7-mL tubes 	–80°C	
Final solution (media)	Oligodendrocyte differentiation	100 ng/mL	Dilute stock solution 1:1000	4°C	



McGill



Bone morphogenetic protein 4 (BMP4)



GenScript #Z02750

• Regulates growth, differentiation, chemotaxis, and apoptosis of various cell types

	\sim				Į	
	MW/Concen	tration	Preparation instruct	ions	Storage	Important information
As supplied	~13.3 k[Da	 Powder; briefly centrifu before opening 	ıge vial	–20°C	 Stable at 4°C but should be kept at -20°C for long-term storage
Stock solution	50 μg/m	۱L	 Dissolve 100 mg in 2 mL sterile ddH₂O cor 0.1% BSA Prepare 100-µL aliquor 0.7-mL tubes 	-	–80°C	 Stable at 4°C for up to 1 week Avoid repeated freeze- thaw cycles
Final solution (media)	Microglia differentiation	50 ng/mL	Dilute stock solution 1:	1000	4°C	





Bone morphogenetic protein 9 (BMP9)/GDF-2



Peprotech #120-07

• Regulates cell differentiation and survival

	~~		Į	
	MW/Concentratio	n Preparation instructions	Storage	Important information
As supplied	24.1 kDa	 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	10 μg/mL	 Dissolve 200 µg in 20 mL sterile ddH₂O containing 0.1% BSA Prepare 100-µL aliquots in 0.7- mL tubes 	–80°C	
Final solution (media)	Cholinergic neuron n differentiation	10 g/mL • Dilute stock solution 1:1000	4°C	







Bovine serum albumin (BSA)



Invitrogen #15260-037

• Supports growth of human hematopoietic progenitor cells in serum-free media formulations

	MW/Concer	ntration	Preparation instructions	Storage	Important information
As supplied/ stock solution	75 mg/mL i	n PBS	• Liquid	–20°C	
Final solution (media)	DA NPC and cortical NPC induction	1 mg/mL	Dilute stock solution 1:75	4°C	







Brain-derived neurotrophic factor (BDNF)



Peprotech #450-02

• Neurotrophic growth factor that supports neuronal proliferation and survival

	\sim			Į	
	MW/Concentra	ition	Preparation instructions	Storage	Important information
As supplied	27.0 kDa		 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	20 µg/mL		 Dissolve 500 µg in 25 mL sterile ddH₂O containing 0.1% BSA Prepare 100-µL aliquots in 0.7- mL tubes 	–80°C	
Final solution	DA and cortical neuron differentiation, forebrain organoid culture	20 ng/mL	Dilute stock solution 1:1000	4°C	
(media)	Motor neuron differentiation, midbrain organoid culture	10 ng/mL	Dilute stock solution 1:2000	4°C	







CD200

Novoprotein #C311

• Immunoregulatory protein



	~~		Į	
	MW/Concentration	Preparation instructions	Storage	Important information
As supplied	30 kDa	 Powder; briefly centrifuge vial before opening 	–20°C	 Stable at room temperature for up to 3 weeks
Stock solution	100 µg/mL	 Dissolve 200 µg in 2 mL sterile ddH₂O containing 0.1% BSA Prepare 100-µL aliquots in 0.7-mL tubes 	–80°C for up to 3 months	 Do not mix by vortex or pipetting Stable at 4°C for 2 to 7 days Minimize freeze-thaw cycles
Final solution (media)	Microglia 100 ng/mL differentiation	Dilute stock solution 1:1000	4°C	







CHIR-99021

Selleckchem #S2924

- Inhibits GSK-3 α/β leading to activation of GSK and insulin signaling

	\sim			Į	
	MW/Concentrat	ion	Preparation instructions	Storage	Important information
As supplied	501.4 g/mol		 Powder; briefly centrifuge vial before opening 	–20°C for up to 3 years	
Stock solution	3 mM		 Dissolve 25 mg in 16.607 mL DMSO Prepare 200-µL aliquots in 1.5-mL tubes 	–80°C for up to 2 years	
Final	DA NPC and cortical NPC induction; motor NPC induction (step 1 and 3 media)	3 µM	Dilute stock solution 1:1000	4°C	
solution (media)	Motor NPC induction (step 2 media)	1 µM	 Dilute stock solution 1:3000 	4°C	
	Midbrain organoid culture	0.8 µM	 Dilute stock solution 1:3750 	4°C	



McGill





Preparation instructions Important information Storage **MW/Concentration** · Powder; briefly centrifuge vial As 22.8 kDa -20°C supplied before opening • Dissolve 250 µg in 25 mL sterile ddH₂O Stock 10 µg/mL containing 0.1% BSA -80°C solution Prepare 100-µL aliquots in 0.7-mL tubes Final Motor neuron 10 Dilute stock solution 1:1000 4°C solution and astrocyte ng/mL (media) differentiation



Peprotech #450-13

• Neural factor that supports neuronal survival

Institut-Hôpital neurologique de Montréal Montreal Neurological institute-Hospital







Compound E

STEMCELL Technologies #73954

Inhibits γ-secretase and Notch processing

	~				
	MW/Concentra	tion	Preparation instructions	Storage	Important information
As supplied	490.5 g/mol	I	 Powder; briefly centrifuge vial before opening 	–20°C for up to 12 months For long-term storage, store with a desiccant	Protect from light
Stock solution	0.1 mM		 Dissolve 1 mg in 20.388 mL DMSO Prepare 100-µL aliquots in 0.7-mL tubes 	–20°C	Avoid repeated freeze- thaw cyclesProtect from light
Final solution (media)	DA neuron, cortical neuron, and motor neuron differentiation	0.1 µM	Dilute stock solution 1:1000	4°C	Protect from light









CX3CL1



Peprotech #300-31

• Chemokine that chemoattracts specific cell types, including microglia cells

	~	,			
	MW/Concer	ntration	Preparation instruction	s Storage	Important information
As supplied	8.5 kD	Da	Powder; briefly centrifuge before opening	vial –20°C	
Stock solution	100 µg/	'nL	 Dissolve 500 μg in 5 mL sterile ddH₂O contair 0.1% BSA Prepare 100-μL aliquots in mL tubes 	-80°C	
Final solution (media)	Microglia differentiation	100 ng/mL	Dilute stock solution 1:100	0 4°C	







DAPT

Selleckchem #S2215

Inhibits γ-secretase



	MW/Concentrat	ion	Preparation instructions	Storage	Important information
As supplied	432.46 g/mol		 Powder; briefly centrifuge vial before opening 	–20°C for up to 3 years	
Stock solution	10 mM		 Dissolve 10 mg in 2.312 mL DMSO Prepare 100-µL aliquots in 0.7-mL tubes 	–80°C for up to 2 years	
Final solution (media)	DA neuron and sensory neuron differentiation	10 µM	Dilute stock solution 1:1000	4°C	







Dexamethasone

Sigma #D4902

Activates glucocorticoid receptor signaling



	\sim				
	MW/Concentrat	tion	Preparation instructions	Storage	Important information
As supplied	392.46 g/mol	I	 Powder; briefly centrifuge vial before opening 	4°C	
Stock solution	10 mM		 Dissolve 25 mg in 6.37 mL DMSO Prepare 20-µL aliquots in 0.7-mL tubes 	-80°C	 Avoid repeated freeze- thaw cycles
Final solution (media)	PBMC reprogramming	1 µM	Dilute stock solution 1:10,000	4°C	





EDDU-001-02

March 2020



Dibutyryl-cAMP (db-cAMP)



Carbosynth #ND07996

• cAMP analog that actives cAMP-dependent kinases

	~				Į	
	MW/Concentra	tion		Preparation instructions	Storage	Important information
As supplied	491.37 g/mo	l	•	Powder; briefly centrifuge vial before opening	–20°C	
Stock solution	0.5 M		•	Dissolve 1 g in 4.07 mL sterile ddH ₂ O Prepare 200-µL aliquots in 0.7-mL tubes	–80°C	
Final solution	DA neuron and cortical neuron differentiation	0.5 mM	•	Dilute stock solution 1:1000	4°C	
(media)	Midbrain organoid culture	125 µM	•	Dilute stock solution 1:4000	4°C	







Dorsomorphin

Tocris #3093

• Inhibits AMPK and type I BMP receptors

	~				
	MW/Concentra	ation	Preparation instructions	Storage	Important information
As supplied	472.4 g/mo	I	 Powder; briefly centrifuge vial before opening 	–20°C	 Product is hygroscopic. Desiccate upon arrival Can be stored at room temperature for up to 6 months
Stock solution	20 mM		 Dissolve 10 mg in 4.237 mL DMSO Prepare 100-µL aliquots in 0.7-mL tubes 	–80°C	 Can be stored at –20°C for up to 1 month
Final	NPC induction	2 µM	Dilute stock solution 1:10,000	4°C	
solution (media)	Forebrain organoid culture	10 µM	Dilute stock solution 1:2000	4°C	





EDDU-001-02

March 2020





Dorsomorphin homolog 1 (DMH1)

Selleckchem #S7146

• Inhibits ALK2, a type 1 BMP receptor

	\sim				
	MW/Concentratio	n Preparatio	on instructions	Storage	Important instructions
As supplied	380.44 g/mol	 Powder; brivial before 	iefly centrifuge opening	–20°C for up to 3 years	
Stock solution	4 mM	 Dissolve 25 16.425 mL Prepare 10 0.7-mL tube 	DMSO 0-µL aliquots in	–80°C for up to 2 years	
Final solution (media)	Motor neuron and astrocyte differentiation	2 µM • Dilute stocł	k solution 1:2000	4°C	





EDDU-001-02



DMEM/F12 media

Gibco #10565-018

• Basal media for supporting growth of mammalian cells

	~~		-		
	MW/Concentration	Preparation	instructions	Storage	Important instructions
As supplied	1x	Solution	1	4°C for up to 2 months from date of manufacture	Protect from light
Final solution (coating)	Laminin solution		tic-Antimycotic as per protocol	4°C	Protect from light
Final solution (media)	iPSC culture; DA neuron, cortical neuron, motor neuron, and astrocyte differentiation; cerebral organoid culture	• Add suppler 1x protocol	nents as per	4°C	 Protect from light







DMEM



Wisent Bioproducts #319-005-CL

• Basal media for supporting growth of mammalian cells

	~~				
	MW/Concent	tration	Preparation instructions	Storage	Important instructions
As supplied	1x		Solution	4°C	
Final solution (media)	MEF culture	1x	 Add supplements as per protocol 	4°C	








Embryonic stem-cell fetal bovine serum (FBS)



Gibco #10439024

• Sustains undifferentiated cellular morphology of embryonic stem cells









Epidermal growth factor (EGF)



Peprotech #AF-100-15

• Growth factor that stimulates growth of epidermal and epithelial cells

	~			L	
	MW/Concentra	tion	Preparation instructions	Storage	Important information
As supplied	6.2 kDa		 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	10 µg/mL		 Dissolve 1 mg in 100 mL sterile ddH₂O containing 0.1% BSA Prepare 100-µL aliquots in 0.7-mL tubes 	–80°C	
Final solution (media)	DA neuron, cortical neuron, and astrocyte differentiation [*]	10 ng/mL	Dilute stock solution 1:1000	4°C	
	Astrocyte differentiation*	20 ng/mL	Dilute stock solution 1:500	4°C	

*Final solution concentration for astrocyte differentiation depends on protocol.







Erythropoietin (EPO)



Peprotech #100-64

• Hormone that stimulates proliferation and differentiation of erythroid progenitor cells

	~~				
	MW/Concentrati	ion	Preparation instructions	Storage	Important information
As supplied	150 U/µg*		 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	2000 U/µL		 Dissolve 13.3 µg in 1 mL 1x PBS Prepare 50-µL aliquots in 0.7-mL tubes 	-80°C	
Final solution (media)	PBMC reprogramming	2 U/mL	Dilute stock solution 1:1000	4°C	

*1 U is equivalent to 6.7 ng.







Essential 8 supplement

Gibco #A1517001 (component #A1517101)

• Supplement for Essential 8 basal medium

	~					
	MW/Concentration	Preparation instructions	Storage	Important instructions		
As supplied	50x	Solution	–20°C for up to 12 months from date of manufacture	 Store in non-frost-free freezer Do not re-freeze thawed solution 		
Final solution (media)	iPSC culture 1x	 Thaw Essential 8 50x supplement for about 1 hour at room temperature or overnight at 4°C. Mix gently Add entire bottle (10 mL) of Essential 8 50x supplement to 500 mL of Essential 8 basal media Prepare 40-mL aliquots in 50- mL conical tubes, seal with Parafilm, and store at -20°C 	4°C for up to 2 weeks –20°C for up to 6 months	 Do not thaw/warm supplement or complete media in a 37°C water bath Once supplement is thawed, use immediately or store at 4°C for up to 2 days Thaw complete media aliquots overnight at 4°C. Do not refreeze aliquots after thawing Warm complete media at room temperature 		









Essential 8 basal media

 \bigcirc

Institut-Hôpital neurologique de Montréal Montreal Neurological Institute-Hospital McGill



Gibco #A1517001 (component #A1516901)

• Xeno-free and feeder-free medium used for the growth and expansion of human iPSCs

	~~			
	MW/Concentration	Preparation instructions	Storage	Important instructions
As supplied	1x	Solution	4°C for up to 12 months from date of manufacture	Protect from light
Final solution (media)	iPSC 1x culture	 Thaw Essential 8 50x supplement for about 1 hour at room temperature or overnight at 4°C. Mix gently Add entire bottle (10 mL) of Essential 8 50x supplement to 500 mL of Essential 8 basal media Prepare 40-mL aliquots in 50- mL conical tubes, seal with Parafilm, and store at -20°C 	4°C for up to 2 weeks –20°C for up to 6 months	 Do not thaw/warm supplement or complete media in a 37°C water bath Once supplement is thawed, use immediately or store at 4°C for up to 2 days Thaw complete media aliquots overnight at 4°C. Do not refreeze aliquots after thawing Warm complete media at room temperature



Fetal bovine serum (FBS)

ThermoFisher #12484028

- Supplement for culture media
- Cryopreservation media



	~~				
	MW/Concentration		Preparation instructions	Storage	Important instructions
As supplied/ stock solution	100%		Solution	–20°C	
Final solution (cryo- preservation)	iPSCs, DA NPCs, cortical NPCs, and motor NPCs	90%	• Add DMSO at 10%	Room temperature	
Final solution (media)	MEF culture	10%	 Dilute stock solution 1:10 	4°C	







Fibroblast growth factor-basic (FGF-b)/FGF-2



Peprotech #100-18B

• Growth factor that promotes proliferation and differentiation of a variety of cell types

	MW/Concentration	on	Preparation instructions	Storage	Important information
As supplied	17.2 kDa		 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	10 µg/mL		 Dissolve 1 mg in 100 mL sterile ddH₂O containing 0.1% BSA Prepare 100-µL aliquots in 0.7-mL tubes 	–80°C	
Final	PBMC reprogramming; DA neuron, cortical neuron, and astrocyte differentiation [*]	10 ng/mL	Dilute stock solution 1:1000	4°C	
solution (media)	Astrocyte differentiation*	20 ng/mL	Dilute stock solution 1:500	4°C	
	Cerebral organoid culture	4 ng/mL	Dilute stock solution 1:2500	4°C	

*Final solution concentration for astrocyte differentiation depends on protocol.







Fibroblast growth factor 8 (FGF-8)



Peprotech #100-25

• Growth factor that promotes cell proliferation and differentiation

	\sim				
	MW/Concentra	tion	Preparation instructions	Storage	Important information
As supplied	22.5 kDa		 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	100 μg/mL		 Dissolve 500 µg in 5 mL sterile ddH₂O containing 0.1% BSA Prepare 100-µL aliquots in 0.7- mL tubes 	80°C	
Final solution (media)	NPC induction and midbrain organoid culture	100 ng/mL	Dilute stock solution 1:1000	4°C	







5-Fluoro-2'-deoxyuridine

Sigma #F0503

· Antineoplastic agent that inhibits thymidylate synthetase



	~~			
	MW/Concentration	Preparation instructions	Storage	Important information
As supplied	246.19 g/mol	 Powder; briefly centrifuge vial before opening 	Room temperature	
Stock solution	100 mM	 Dissolve 100 mg in 406 µL sterile ddH₂O Prepare 20-µL aliquots in 0.7-mL tubes 	–80°C	 Stable at 4°C for up to 2 weeks Protect from light
Final solution (media)	DA neuron, cortical neuron, and motor 1 μ neuron differentiation	M • Dilute stock solution 1:100,000	4°C	Protect from light







Gentle Cell Dissociation Reagent



STEMCELL Technologies #07174

• Enzyme-free reagent for dissociation of human embryonic stem cells or human (iPSCs) cells into cell aggregates for routine passaging or into a single-cell suspension









Glial-derived neurotrophic factor (GDNF)



Peprotech #450-10

• Neurotrophic factor that promotes dopamine uptake and survival and morphological differentiation of midbrain neurons

	~~				
	MW/Concentra	ition	Preparation instructions	Storage	Important information
As supplied	30.4 kDa		 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	20 µg/mL		 Dissolve 500 µg in 25 mL sterile ddH₂O containing 0.1% BSA Prepare 100-µL aliquots in 0.7- mL tubes 	–80°C	
Final	DA neuron and cortical neuron differentiation	20 ng/mL	Dilute stock solution 1:1000	4°C	
solution (media)	Midbrain organoid culture	10 ng/mL	Dilute stock solution 1:2000	4°C	




McGill



GlutaMAX-I

Gibco #35050-061

 Supplement for adherent and suspension culture of mammalian cells that is an alternative to L-glutamine, with increased stability that improves cell health













Heparin

Sigma #H3149

• Enhances the antithrombin-mediated inactivation of proteases in the coagulation pathway

	~~						
	MW/Concentration	on		Preparation instructions	St	torage	Important information
As supplied	Mixture of polyanion chai 6–30 kDa	ins of MW	•	Powder; briefly centrifuge vial before opening	_	-20°C	
Stock solution	20 mg/mL		•	Dissolve 477 mg in 23.85 mL 1x PBS Prepare 50-µL aliquots in 0.7-mL tubes	-	-80°C	
Final solution	Astrocyte differentiation	2 µg/mL	•	Dilute stock solution 1:10,000		4°C	
(media)	Midbrain and cerebral organoid culture	1 µg/mL	•	Dilute stock solution 1:20,000		4°C	









Culture Reagents: I–P





EDDU-001-02

March 2020



Insulin

Sigma #I2643

• Regulates cellular uptake, utilization, and storage of glucose, amino acids, and fatty acids and inhibits breakdown of glycogen, protein, and fat

	~~				
	MW/Concentra	tion	Preparation instructions	Storage	Important information
As supplied	5807.57 g/m	ol	 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	5 mg/mL		 Dissolve 50 mg in 10 mL sterile ddH₂O containing 0.1% BSA Prepare 200-µL aliquots in 0.7-mL tubes 	–80°C	
Final solution	Midbrain organoid culture	2.5 µg/mL	Dilute stock solution 1:2000	4°C	
(media)	Cerebral organoid culture	1.25 µg/mL	Dilute stock solution 1:4000	4°C	









Insulin-like growth factor 1 (IGF-1)



Peprotech #100-11

• Growth factor that stimulates proliferation and differentiation of various cell types

	~~			l	
	MW/Concentra	tion	Preparation instructions	Storage	Important information
As supplied	7.6 kDa		 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	10 µg/mL		 Dissolve 1 mg in 100 mL sterile ddH₂O containing 0.1% BSA Prepare 100-µL aliquots in 0.7- mL tubes 	–80°C	
Final solution (media)	PBMC reprogramming, motor neuron differentiation	10 ng/mL	Dilute stock solution 1:1000	4°C	








Insulin-transferrin-sodium selenite supplement (ITS) Gibco #41400045

• Insulin promotes glucose and amino acid uptake, lipogenesis, intracellular transport, and protein and nucleic acid synthesis. Transferrin is an iron carrier and may reduce toxic levels of oxygen radicals and peroxide. Selenium, as sodium selenite, is a cofactor for glutathione peroxidase and other proteins, and an antioxidant












Interleukin-3 (IL-3)

Peprotech #200-03

• Hematopoietic growth factor that promotes the survival, differentiation and proliferation of specific types of committed progenitor cells, including those of the erythroid lineage





McGill



Interleukin-6 (IL-6)

Peprotech #200-06

Cytokine with diverse biological functions



	~~			
	MW/Concentration	Preparation instructions	Storage	Important information
As supplied	20.9 kDa	 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	50 µg/mL	 Dissolve 500 μg in 10 mL sterile ddH₂O containing 0.1% BSA Prepare 100-μL aliquots in 0.7-mL tubes 	–80°C	
Final solution (media)	0	50 • Dilute stock solution 1:1000 /mL	4°C	







Interleukin-34 (IL-34)

Peprotech #200-34

• Ligand for colony-stimulating factor-1 receptor (CSF1R)



	~	/			
	MW/Conce	ntration	Preparation instructions	Storage	Important information
As supplied	52.5 k	Da	 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	100 µg/mL		 Dissolve 500 µg in 5 mL sterile ddH₂O containing 0.1% BSA Prepare 100-µL aliquots in 0.7-mL tubes 	–80°C	
Final solution (media)	Microglia differentiation	100 ng/mL	Dilute stock solution 1:1000	4°C	









KnockOut serum replacement



Gibco #10828028

Supplement for KnockOut DMEM/F12 basal media

	\sim				1
	MW/Concentrat	ion	Preparation instructions	Storage	Important information
As supplied	100%		 Solution Prepare 10-mL aliquots in 15-mL tubes 	–20°C for up to 18 months from date of manufacture	 Avoid repeated freeze- thaw cycles Can be stored at 4°C for up to 4 weeks Protect from light
Final solution (media)	PBMC reprogramming, iPSC culture, cerebral organoid culture	20%	 Thaw aliquot of KnockOut serum replacement overnight at 4°C Add 20 mL KnockOut serum replacement to 100 mL KnockOut DMEM/F12 Add supplements as per protocol 	4°C for up to 10 days	 Do not thaw/warm serum replacement or complete media in a 37°C water bath Warm only the volume of complete media required for that day's use Protect complete media from light







KnockOut DMEM/F12

Gibco #12660012

 Low osmolality medium without L-glutamine or HEPES buffer for the growth of human embryonic stem cells and iPSCs

	MW/Concentration		Preparation instructions	Storage	Important information
As supplied	1x		Solution	4°C	
Final solution (media)	PBMC reprogramming, iPSC culture	1x	 Thaw aliquot of KnockOut serum replacement overnight at 4°C Add 20 mL KnockOut serum replacement to 100 mL KnockOut DMEM/F12 Add supplements as per protocol 	4°C for up to 10 days	 Do not thaw/warm serum replacement or complete media in a 37°C water bath Warm only the volume of complete media required for that day's use Protect complete media from light







Laminin



Invitrogen #23017-015

• Extracellular matrix protein that supports adhesion, proliferation, and differentiation of many cell types

	~~				
	MW/Concentration	n	Preparation instructions	Storage	Important information
As supplied/ stock solution	0.5–2.0 mg/mL [*] in 50 mM (pH 7.4), 0.15 M NaCl	Tris-HCI	 Solution Thaw solution slowly at 4°C Prepare 100-µL aliquots in 0.7-mL polypropylene microcentrifuge tubes 	–80°C	 If frozen solution is warmed too quickly, laminin may form a gel and cannot be reactivated for use Do not freeze-thaw repeatedly
Final solution (coating)	DA neuron, cortical neuron, motor neuron, and astrocyte differentiation	5 μg/mL	 Dilute stock solution to 5 µg/mL in DMEM/F12 	4°C for up to 2 weeks	 Do not store for more than 2 weeks
Final solution (media)	DA neuron, cortical neuron, and motor neuron differentiation	1 µg/mL	 Dilute stock solution to 1 μg/mL 	4°C for up to 2 weeks	 Do not store for more than 2 weeks

*Verify exact concentration on tube label.









Laminin

Sigma #L2020

• Extracellular matrix protein that supports adhesion, proliferation, and differentiation of many cell types

-	~					
	MW/Concentration			Preparation instructions	Storage	Important information
As supplied/ stock solution	~1 mg/mL [*] in 50 mM Tris-HCl (pH 7.5 150 mM NaCl	5),	•	Solution Thaw solution slowly at 4°C Prepare 100-µL aliquots in 0.7-mL polypropylene microcentrifuge tubes	–80°C	 If the frozen solution is warmed too quickly, laminin may form a gel and cannot be reactivated for use
Final solution (coating)	DA NPC and cortical NPC induction; DA neuron, cortical neuron, motor neuron, and astrocyte differentiation	5 µg/mL	•	Dilute stock solution to 5 μg/mL in DMEM/F12	4°C for up to 2 weeks	 Do not use if discoloration or web formations appear on coated surfaces
Final	DA neuron, cortical neuron, and motor neuron differentiation	1 µg/mL	•	Dilute stock solution to 1 μg/mL	4°C for up to 2 weeks	Do not store for more than 2 weeks
solution (media)	Midbrain organoid culture	200 ng/mL	•	Dilute stock solution to 200 ng/mL	4°C for up to 2 weeks	 Do not store for more than 2 weeks

*Verify exact concentration on tube label.







L-Ascorbic acid 2-phosphate (AA2P; sesquimagnesium salt hydrate) Sigma #A8960



• Long-acting ascorbic acid derivative that stimulates collagen expression and formation

	\sim				
	MW/Concent	ration	Preparation instructions	Storage	Important information
As supplied	289.54 g/mol		 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	32 mg/mL		 Dissolve 32 mg in 1 mL sterile ddH₂O Prepare 100-µL aliquots in 0.7-mL tubes 	4°C	
Final solution (media)	PBMC reprogramming	64 μg/mL	Dilute stock solution 1:500	4°C	







LDN193189

Sigma #SML0559

• Dorsomorphin derivative that inhibits ALK2 and ALK3



	\sim			
	MW/Concentration	Preparation instructions	Storage	Important information
As supplied	406.48 g/mol	 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	1 mM	 Dissolve 5 mg in 12.3 mL sterile ddH₂O containing 0.1% BSA Prepare 100-µL aliquots in 0.7-mL tubes 	–80°C	
Final solution (media)	Astrocyte 0.1 μ differentiation	• Dilute stock solution 1:10,000	4°C	







L-Glutamine

Wisent Bioproducts #609-065-EL

 Amino acid required for cell culture involved in the formation of purine and pyrimidine nucleotides, amino sugars, glutathione, L-glutamate, other amino acids; protein synthesis; and glucose production

	\sim				
	MW/Concentrat	ion	Preparation instructions	Storage	Important information
As supplied/ stock solution	200 mM		Solution	–20°C	
Final solution (media)	PBMC reprogramming, MEF culture	2 mM	 Dilute stock solution 1:100 	4°C	







Lithium chloride (LiCl)

Sigma #L7026

• Greatly enhances the generation of iPSCs

	~~		Į	
	MW/Concentrati	on Preparation instructions	Storage	Important information
As supplied	8 M in H ₂ O	Solution	4°C	
Stock solution	2 M	 Dilute in sterile ddH₂O to achieve 2 M Filter solution Prepare 100-µL aliquots in 0.7- mL tubes 	4°C	
Final solution (media)	PBMC reprogramming	2 mM • Dilute stock solution 1:1000	4°C	








Macrophage colony stimulating factor (M-CSF)



Peprotech #300-25

• Hematopoietic factor that regulates cellular proliferation, differentiation, and survival in blood monocytes, tissue macrophages, and their respective progenitor cells










Matrigel growth factor reduced (GFR) basement membrane matrix

Corning #356230

 Solubilized basement membrane extract containing laminin (a major component), collagen IV, heparin sulfate proteoglycans, entactin/nidogen, and a number of growth factors










Matrigel human embryonic stem cell (hESC)– qualified matrix Corning Millipore #354277

 Soluble basement membrane extract containing laminin, collagen IV, entactin and heparan sulfate, and proteoglycan



*Exact concentration on tube label.








2-Mercaptoethanol (β-mercaptoethanol [BME])



Gibco #21985023

• Reducing agent used in culture media to prevent toxic levels of oxygen radicals







2-Mercaptoethanol (β-mercaptoethanol [BME])



Merck #8057400005

• Reducing agent used in culture media to prevent toxic levels of oxygen radicals

	~~	•			
	MW/Concen	tration	Preparation instructions	Storage	Important information
As supplied/ stock solution	78.12 g/mol (14	4.25 mM)	Solution	Room temperature	
Final	Midbrain organoid culture	0.00035%	 Dilute stock solution to 0.00035% 	4°C	
solution (media)	Cerebral organoid culture	0.00035% or 0.0007%	 Dilute stock solution to 0.00035% or 0.0007%, as per protocol 	4°C	







MEM nonessential amino acid (NEAA) solution



Wisent #321-011-EL

• Supplement to increase cell growth and viability

	~~					
	MW/Concentration			Preparation instructions	Storage	Important information
As supplied/ stock solution	100x		•	Solution	4°C	
Final solution (media)	PBMC reprogramming; DA NPC and cortical NPC induction; DA neuron, cortical neuron, and astrocyte differentiation; midbrain and cerebral organoid culture [*]	1x	•	Dilute stock solution 1:100	4°C	
	Cerebral organoid culture [*]	0.5x	•	Dilute stock solution 1:200	4°C	

*Final solution concentration for cerebral organoid culture depends on cerebral organoid stage, as per protocol.







Mitomycin C



Sigma #M4287

• Antibiotic and double-stranded DNA alkylating agent that inhibits cell proliferation

	\sim		Į	
	MW/Concentratio	on Preparation instructions	Storage	Important information
As supplied	334.33 g/mol	 Powder; briefly centrifuge vial before opening 	4°C	
Stock solution	1 mg/mL	 Dissolve 2 mg in 2 mL sterile ddH₂O Prepare 50-µL aliquots in 0.7- mL tubes 	–80°C	 Can be stored at 4°C for up to 1 week (pH 6–9) If a precipitate forms, prepare a fresh stock solution Protect from light
Final solution (media)	DA neuron and cortical neuron differentiation	1 μg/mL • Dilute stock solution 1:1000	4°C	Protect from light







mTeSR1 5x supplement



STEMCELL Technologies #85850 (component #85852)

• Supplement for mTeSR1 basal media

	\sim		Į	
	MW/Concentration	Preparation instructions	Storage	Important instructions
As supplied	5x	Solution	–20°C for up to 3 months	
Final solution (media)	iPSC culture 1x	 Thaw mTeSR1 5x supplement overnight at 4°C. Mix thoroughly Add 100 mL mTeSR1 5x supplement to 400 mL mTeSR1 basal media to obtain 500 mL of complete media. Mix well Prepare 40-mL aliquots in 50- mL conical tubes, seal with Parafilm, and store at -20°C 	4°C for up to 2 weeks –20°C for up to 6 months	 Do not thaw/warm supplement or complete media in a 37°C water bath Once supplement is thawed, use immediately or store at 4°C for up to 2 days Thaw complete media aliquots overnight at 4°C. Do not refreeze aliquots after thawing Warm complete media at room temperature







mTeSR1 basal media



STEMCELL Technologies #85850 (component #85851)

• Feeder-free cell culture media for human embryonic stem cells and iPSCs

As	MW/Concentra	luon	Preparation instructions Solution	Storage 4°C	Important instructions	
supplied	1x		Thaw mTeSR1 5x supplement	4 0	• Do not thaw/warm supplement	
Final solution (media)	iPSC culture	1x	 overnight at 4°C. Mix thoroughly Add 100 mL mTeSR1 5x supplement to 400 mL mTeSR1 basal media to obtain 500 mL of complete media. Mix well Prepare 40-mL aliquots in 50- mL conical tubes, seal with Parafilm, and store at -20°C 	4°C for up to 2 weeks –20°C for up to 6 months	 or complete media in a 37°C water bath Once supplement is thawed, use immediately or store at 4°C for up to 2 days Thaw complete media aliquots overnight at 4°C. Do not refreeze aliquots after thawing Warm complete media at room temperature 	







N-2 supplement

Gibco #17502048

• Supplement that promotes neuronal cell growth



	~~				
	MW/Concentration		Preparation instructions	Storage	Important information
As supplied/ stock solution	100x		Solution	–20°C	Protect from light
Final solution (media)	DA neuron, cortical neuron, and astrocyte differentiation; midbrain organoid culture	1x	 Dilute stock solution 1:100 	4°C	Protect from light
	Motor neuron differentiation, cerebral organoid culture	0.5x	 Dilute stock solution 1:200 	4°C	Protect from light







β-Nerve growth factor (β-NGF)

Peprotech #450-01

• Neurotrophic factor that is crucial for development and preservation of sensory and sympathetic nervous systems

MW/Concentration		Preparation instructions	E Storage	Important information
13.5 kDa		 Powder; briefly centrifuge vial before opening 	–20°C	
100 μg/mL		 Dissolve 500 µg in 5 mL sterile ddH₂O containing 0.1% BSA Prepare 100-µL aliquots in 0.7-mL tubes 	–80°C	
Cholinergic neuron differentiation	100 ng/mL	Dilute stock solution 1:1000	4°C	
	13.5 kDa 100 μg/mL Cholinergic neuron	13.5 kDa 100 μg/mL Cholinergic neuron 100	 13.5 kDa Powder; briefly centrifuge vial before opening Dissolve 500 μg in 5 mL sterile ddH₂O containing 0.1% BSA Prepare 100-μL aliquots in 0.7-mL tubes Cholinergic neuron 100 	13.5 kDa • Powder; briefly centrifuge vial before opening -20°C 100 μg/mL • Dissolve 500 μg in 5 mL sterile ddH ₂ O containing 0.1% BSA -80°C • Prepare 100-μL aliquots in 0.7-mL tubes -80°C








Neurobasal (NB) media



Life Technologies #21103-049

Basal media for long-term maintenance and maturation of pure pre-natal and embryonic neuronal cell populations









Neurotrophin-3 (NT3)



Peprotech #450-03

• Neurotrophic factor that promotes the growth and survival of nerve and glial cells

	MW/Concentr	ation		Preparation instructions	Storage	Important information
As supplied	13.6 kDa		•	Powder; briefly centrifuge vial before opening	–20°C	
Stock solution	10 µg/mL		•	Dissolve 50 µg in 5 mL sterile ddH ₂ O containing 0.1% BSA Prepare 100-µL aliquots in 0.7-mL tubes	-80°C	
Final solution	Sensory neuron and oligodendrocyte differentiation	10 ng/mL	•	Dilute stock solution 1:1000	4°C	
(media)	Forebrain organoid 20 ng/mL culture		•	Dilute stock solution 1:500	4°C	

EDDU-001-02

March 2020



McGill

Noggin

Peprotech #120-10C

Inhibits TGF-β ligands during developmental processes



	~~					
	MW/Concentrati	on		Preparation instructions	Storage	Important information
As supplied	46 kDa		•	Powder; briefly centrifuge vial before opening	–20°C	
Stock solution	200 μg/mL		•	Dissolve 500 µg in 2.5 mL sterile ddH ₂ O containing 0.1% BSA Prepare 100-µL aliquots in 0.7-mL tubes	–80°C	
Final solution (media)	NPC induction, midbrain organoid culture	200 ng/mL	•	Dilute stock solution 1:1000	4°C	







Penicillin-Streptomycin

Wisent #450-200-EL

• Contains the antibiotics penicillin and streptomycin



	~~		Į	
	MW/Concentration	Preparation instructions	Storage	Important information
As supplied/ stock solution	100X	Solution	–20°C	
Final solution (media)	PBMC reprogramming, MEF culture, midbrain and cerebral organoid culture	• Dilute stock solution 1:100	4°C	







Platelet-derived growth factor-AA (PDGF-AA)



Peprotech #100-13A

Mitogen involved in a number of biological processes, including embryonic neuron development

	~				
	MW/Concentra	tion	Preparation instructions	Storage	Important information
As supplied	28.5 kDa		 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	10 µg/mL		 Dissolve 50 µg in 5 mL sterile ddH₂O containing 0.1% BSA Prepare 100-µL aliquots in 0.7- mL tubes 	-80°C	
Final solution (media)	Oligodendrocyte differentiation	10 ng/mL	Dilute stock solution 1:1000	4°C	







Phosphate-buffered saline (PBS)

Wisent #311-010-CL

• Balanced salt solution used for a variety of cell culture applications, such as washing cells and preparing reagents













2-Phospho-L-ascorbic acid (AA2P) trisodium salt



Sigma #49752

• Long-acting ascorbic acid derivative that stimulates collagen expression and formation

	\sim				
	MW/Concentra	tion	Preparation instructions	Storage	Important information
As supplied	322.05 g/mo	I	 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	50 mg/mL		 Dissolve 500 mg in 10 mL sterile ddH₂O Filter solution Prepare 100-µL aliquots in 0.7- mL tubes 	-80°C	
Final solution (media)	PBMC reprogramming	50 μg/mL	Dilute stock solution 1:1000	4°C	







Poly-L-ornithine (PO)

Sigma #P3655

• Synthetic amino acid polymer coating that acts as a charge modifier

	~~			Į	1
	MW/Concentra	ation	Preparation instructions	Storage	Important information
As supplied	30–70 kDa	l	• Powder	–20°C	
Stock solution	1 mg/mL		 Dissolve 100 mg in 100 mL sterile 1x PBS Prepare 10-mL aliquots in 15-mL tubes 	–20°C	
Final solution (coating)	DA neuron, cortical neuron, motor neuron, and astrocyte differentiation	10 µg/mL	 Dilute stock solution 1:100 in 1x PBS 	4°C for up to 2 months	 Ensure transfer of all stock solution by rinsing tube with 1x PBS twice









Polyvinyl alcohol (PVA)

Sigma #P8136

• Hydrophilic linear polymer which forms copolymers of vinyl alcohol and vinyl acetate

	MW/Concentration		Preparation instructions	Storage	Important information
As supplied	28.5 kDa		• Powder	–20°C	
Stock solution	10 mg/mL		 Dissolve 100 mg in 10 mL ddH₂O Filter solution Prepare 100-µL aliquots in 0.7-mL tubes 	–80°C	
Final solution (coating)	Microglia 10 differentiation	0 μg/mL	Dilute stock solution 1:1000	4°C	







Purmorphamine

Sigma #SML-0868

• Binds to the Smo receptor and activates Hedgehog signaling

	~~				
	MW/Concentratio	on	Preparation instructions	Storage	Important information
As supplied	520.6 g/mol		 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	2 mM		 Dissolve 5 mg in 4.802 mL DMSO Prepare 100-µL aliquots in 0.7-mL tubes 	–80°C	 Prepare as accurately as possible because the working concentration range is very narrow
	DA neuron differentiation	2 µM	Dilute stock solution 1:1000	4°C	
Final solution (media)	Motor NPC induction (step 2 and 3 media)	0.5 µM	Dilute stock solution 1:4000	4°C	 Use the smallest tip and a well-calibrated pipette
	Motor neuron differentiation	0.1 µM	 Dilute stock solution 1:20,000 	4°C	











Culture Reagents: Q–Z





EDDU-001-02

March 2020



Retinoic acid (RA)

Sigma #R2625

• Activates transcription factors that regulate cell growth and differentiation

	~~			Į	
	MW/Concentra	ation	Preparation instructions	Storage	Important information
As supplied	300.44 g/ma	ol	 Powder; briefly centrifuge vial before opening 	–20°C for up to 1 year	Stable if unopened in ampuleProtect from light and air
Stock solution	1 mM		 Dissolve 300 mg in 10 mL DMSO Dilute solution 1:10 in EtOH Prepare 20-µL aliquots in 0.7-mL tubes 	–80°C	 Use all of the powder immediately after opening the ampule Protect from light and air
Final solution (media)	Motor neuron differentiation	0.5 µM	Dilute stock solution 1:2000	4°C	 Always use a new aliquot of stock solution to prepare the
	Motor NPC induction (step 2 and 3 media)	0.1 µM	 Dilute stock solution 1:10,000 	4°C	final solution Protect from light and air









SB431542

Selleckchem #S1067

• Inhibits ALK5 and the TGF- β /Activin/NODAL pathway

	~~				
	MW/Concentrat	ion	Preparation instructions	Storage	Important information
As supplied	384.4 g/mol		 Powder; briefly centrifuge vial before opening 	–20°C for up to 3 years	
Stock solution	10 mM		 Dissolve 10 mg in 2.604 mL DMSO Prepare 100-µL aliquots in 0.7- mL tubes 	–80°C for up to 2 years	
Final solution (media)	DA NPC and cortical NPC induction, astrocyte differentiation [*] , midbrain organoid culture	10 µM	Dilute stock solution 1:1000	4°C	
(media)	Motor neuron and astrocyte 2 µM differentiation [*]		Dilute stock solution 1:5000	4°C	

*Final solution concentration for astrocyte differentiation depends on protocol.







Smoothened Agonist (SAG)

Millipore #566660

• Activates Smo signalling, Smo internalization, and Hedgehog signaling

	\sim)		Į	
	MW/Concent	ration	Preparation instructions	Storage	Important information
As supplied	599 g/mc	bl	 Powder; briefly centrifuge vial before opening 	–20°C	Can be stored at 4°CProtect from light
Stock solution	1 mM		 Dissolve 1 mg in 1.669 mL sterile ddH₂O containing 0.1% BSA Prepare 100-µL aliquots in 0.7-mL tubes 	–80°C	 Can be stored at -20°C for up to 6 months Protect from light
Final solution (media)	NPC induction	0.25 µM	Dilute stock solution 1:4000	4°C	Protect from light







Sodium butyrate

Sigma #B5887

Inhibits histone deacetylases



	\sim)			
	MW/Concentr	ration	Preparation instructions	Storage	Important information
As supplied	110 g/mo	l	 Powder; briefly centrifuge vial before opening 	Room temperature	
Stock solution	125 mM		 Dissolve 250 mg in 18.18 mL sterile ddH₂O Prepare 100-µL aliquots in 0.7-mL tubes 	–80°C	
Final solution (media)	PBMC reprogramming	250 µM	Dilute stock solution 1:500	4°C	





EDDU-001-02





Sonic hedgehog (SHH)



Peprotech #100-45

 Morphogen that activates Hedgehog signaling and is involved in central nervous system patterning during development

	~~			I	
	MW/Concentra	ation	Preparation instructions	Storage	Important information
As supplied	20.0 kDa		 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	200 µg/mL		 Dissolve 500 µg in 2.5 mL sterile ddH₂O containing 0.1% BSA Prepare 100-µL aliquots in 0.7-mL tubes 	–80°C	
Final solution (media)	NPC induction	200 ng/mL	Dilute stock solution 1:1000	4°C	
	Midbrain organoid culture	100 ng/mL	Dilute stock solution 1:2000	4°C	





Sonic hedgehog (SHH; C24II)

GenScript #Z03067

 Morphogen that activates Hedgehog signaling and is involved in central nervous system patterning during development

	~	•				
	MW/Concen	tration	Preparation inst	tructions	Storage	Important information
As supplied	20.0 kE	Da	 Powder; briefly cer before opening 	ntrifuge vial	–80°C for up to 6 months from date of receipt	
Stock solution	200 µg/mL		 Dissolve 500 µg in sterile ddH₂O cont BSA Prepare 100-µL al mL tubes 	aining 0.1%	80°C	 Can be stored at 4°C for up to 2 weeks or at -20°C for up to 3 months
Final solution (media)	NPC induction	200 ng/mL	Dilute stock solution	on 1:1000	4°C	



McGill



Stem cell factor (SCF)



Peprotech #300-07

 Hematopoietic growth factor that is essential for the survival, proliferation and differentiation of hematopoietic cells committed to the melanocyte and germ cell lineages

			l	
	MW/Concentration	Preparation instructions	Storage	Important information
As supplied	18.4 kDa	 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	10 µg/mL	 Dissolve 50 µg in 5 mL sterile ddH₂O containing 0.1% BSA Prepare 100-µL aliquots in 0.7-mL tubes 	–80°C	
Final solution (media)	PBMC 100 ng/mL reprogramming	Dilute stock solution 1:100	4°C	







Stemline II hematopoietic stem cell expansion medium Sigma #S0192



• Hematopoietic stem cell expansion media for differentiated and undifferentiated cells

	~~				
	MW/Concentration	Preparation instructions	Storage	Important information	
As supplied	1x	Solution	4°C	Protect from light	
Final solution (media)	1x	 Add supplements as per protocol 	4°C	Protect from light	








StemPro Accutase Cell Dissociation Reagent

ThermoFisher #A1110501

• Cell detachment solution of proteolytic and collagenolytic enzymes











Thiazovivin

Tocris #3845

Inhibits ROCK and the RHO/ROCK pathway



	~			l	
	MW/Concentratio	n	Preparation instructions	Storage	Important information
As supplied	311 g/mol		 Powder; briefly centrifuge vial before opening 	–20°C for up to 6 months	
Stock solution	2 mM		 Dissolve 10 mg in 16.367 mL DMSO Prepare 100-µL aliquots in 0.7-mL tubes 	–20°C for up to 1 month	
Final solution (media)	iPSC culture; DA neuron, cortical neuron, and motor neuron differentiation	2 µM	Dilute stock solution 1:1000	4°C	







1-Thioglycerol

Sigma #M1753

• Stimulates cell proliferation



	\sim				
	MW/Concentrat	tion	Preparation instructions	Storage	Important information
As supplied	1.25 g/mL (11.57 M)		Solution	4°C	
Stock solution	0.4 M		 Dilute in sterile ddH₂O to achieve 0.4 M Prepare 100-µL aliquots in 0.7-mL tubes 	4°C	
Final solution (media)	PBMC reprogramming	0.4 mM	Dilute stock solution 1:1000	4°C	









Thrombopoietin (TPO)

Peprotech #300-18

• Growth factor that stimulates the proliferation and maturation of megakaryocytes, and promotes increased circulating levels of platelets

	~	,			
	MW/Concer	ntration	Preparation instructions	Storage	Important information
As supplied	18.4 kDa		 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	50 μg/mL		 Dissolve 250 µg in 5 mL sterile ddH₂O containing 0.1% BSA Prepare 100-µL aliquots in 0.7- mL tubes 	–80°C	
Final solution (media)	Microglia differentiation	50 ng/mL	Dilute stock solution 1:1000	4°C	









Transforming growth factor β3 (TGF-β3)

Peprotech #100-36E

• Cytokine that regulates cell proliferation, growth, differentiation and motility, as well as synthesis and deposition of the extracellular matrix

	~~			
	MW/Concentration	Preparation instructions	Storage	Important information
As supplied	25.0 kDa	 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	1 μg/mL	 Dissolve 50 µg in 50 mL sterile ddH₂O containing 0.1% BSA Prepare 100-µL aliquots in 0.7- mL tubes 	-80°C	
Final solution (media)	DA neuron and cortical neuron 1 ng/mL differentiation	Dilute stock solution 1:1000	4°C	









3,3',5-Triiodo-L-thyronine (T3)



Sigma #T6397

• Thyroid hormone that regulates cell differentiation and protein expression

	\sim			
	MW/Concentration	Preparation instructions	Storage	Important information
As supplied	672.96 g/mol	 Powder; briefly centrifuge vial before opening 	–20°C	
Stock solution	60 μg/mL	 Dissolve 6 mg in 10 mL 0.1N NaOH Add 1 mL of solution to 9 mL 1x PBS Filter solution Prepare 100-µL aliquots in 0.7-mL tubes 	–80°C	 Can be stored at 4°C for up to 1 month or at -20°C Avoid repeated freezing and thawing of aliquots
Final solution (media)	Oligodendrocyte 60 differentiation ng/n	Dilute stock solution 1:1000	4°C	







Valproic acid (VPA)

Sigma #P4543

Anticonvulsant; inhibits histone deacetylates and enzymes involved in GABA metabolism

	~				
	MW/Concentr	ation	Preparation instructions	Storage	Important information
As supplied	166.2 g/mo	bl	Powder	Room temperature	
Stock solution	0.5 M		 Dissolve 831 mg in 10 mL ddH₂O Filter solution Prepare 100-µL aliquots in 0.7 mL tubes 	–80°C	
Final solution	DA neuron, cortical neuron, and motor neuron differentiation	0.5 mM	Dilute stock solution 1:1000	4°C	
(media)	Astrocyte differentiation	10 µM	Dilute stock solution 1:50,000	4°C	







XAV939

Sigma #X3004

• Inhibits Tankyrase and Wnt/β-catenin signaling

	~~			
	MW/Concentration	Preparation instructions	Storage	Important information
As supplied	312 g/mol	 Powder; briefly centrifuge vial before opening 	–20°C	Can be stored at room temperature
Stock solution	10 mM	 Dissolve 25 mg in 8.012 mL DMSO Prepare 100-µL aliquots in 0.7-mL tubes 	–80°C	
Final solution (media)	GABAergic neuron 1 differentiation	• Dilute stock solution 1:1000	4°C	







Y-27632

Selleckchem #S1049

• Inhibits ROCK

	~~			I	
	MW/Concentrat	ion	Preparation instructions	Storage	Important information
As supplied	320.26 g/mol		 Powder; briefly centrifuge vial before opening 	–20°C for up to 3 years	
Stock solution	10 mM		 Dissolve 50 mg in 15.612 mL sterile ddH₂O Prepare 200-µL aliquots in 0.7- mL tubes 	-80°C for up to 2 years	
Final solution (media)	iPSC culture; DA neuron, cortical neuron, and motor neuron differentiation; midbrain organoid culture	10 µM	Dilute stock solution 1:1000	4°C	
	Cerebral organoid culture	50 µM	Dilute stock solution 1:200	4°C	











You have reached the end of the protocol.





EDDU-001-02

March 2020

