



Institut-Hôpital  
neurologique de Montréal  
Montreal Neurological  
Institute-Hospital



**EDDU Protocols**

# **Culture Reagents**

Authors: Carol Xiuqing Chen, Chanshuai Han, Meghna Mathur,  
Nguyen-Vi Mohamed, and Cecilia Rocha  
v2.0



# Antibiotic-Antimycotic

## Gibco #15240-062

- Contains the antibiotics penicillin and streptomycin and the antifungal amphotericin B



MW/Concentration

Preparation instructions

As supplied/  
stock  
solution

100x

- Solution

date of manufacture

Final  
solution  
(media)

iPSC culture;  
DA neuron, cortical  
neuron, and motor  
neuron differentiation

1x

- Dilute 1:100

4°C

Final  
solution  
(coating)

Laminin solution

1x

- Dilute 1:100

4°C


Select to go to the **Table of Contents** screen

Select to go to the **Materials** section

Select to return to this **Help** screen



Select to navigate to the **previous screen**  
or to the **next screen**



# Table of Contents



## Introduction

Objectives

Technical and safety considerations

## Materials

Culture reagents: A–H

Culture reagents: I–P

Culture reagents: Q–Z

## Reagents

Culture reagents: A–H

Culture reagents: I–P

Culture reagents: Q–Z





# Introduction



# Introduction



## Objectives

- This document provides key information for reagents used in culture media for:
  - Peripheral blood mononuclear cell (PBMC) reprogramming
  - Induced pluripotent stem cell (iPSC) culture
  - Mouse embryonic fibroblast (MEF) culture
  - Induction of dopaminergic (DA) neural progenitor cells (NPCs), cortical NPCs, and motor NPCs
  - NPC culture
  - Differentiation of dopaminergic (DA) neurons, cortical neurons, motor neurons, cholinergic neurons, sensory neurons, GABAergic neurons, astrocytes, microglia, and oligodendrocytes
  - 3D culture of forebrain, midbrain, and cerebral organoids
- The following information is provided for each reagent as supplied and for stock and final solutions, where applicable:
  - Molecular weight (MW) of solid powders or concentration of solutions
  - Preparation instructions
  - Storage conditions
  - Important information for handling and usage



# Introduction



## Technical and Safety Information

- This information in this document is not intended to replace the product information or material safety data sheet (MSDS) from the manufacturer. Refer to all documentation provided by the manufacturer prior to using any reagent.
- Follow all laboratory safety requirements as outlined by McGill University Environmental Health and Safety (EHS).
- All reagents in this document are for use in cell culture and should be kept sterile. All solutions should be prepared using sterile technique.
- For reagents supplied as solid powders, reconstitute the entire vial or bottle to prepare stock solutions when appropriate.
- Final solution information is provided for the most common application(s) and is subject to change depending on the application and in the event of protocol updates. Prepare final solutions as per the most recent version of the protocol.





# Materials





## Select a culture reagent to go to its screen. Culture Reagents A–H

- [Activin-A](#)
- [Antibiotic-Antimycotic](#)
- [Ascorbic acid \(AA\)](#)
- [B-27 supplement](#)
- [B-27 supplement, minus Vitamin A](#)
- [Biotin](#)
- [Bone morphogenetic protein 4 \(BMP4\)](#)
- [Bone morphogenetic protein 9 \(BMP9\)/GDF-2](#)
- [Bovine serum albumin \(BSA\)](#)
- [Brain-derived neurotrophic factor \(BDNF\)](#)
- [CD200](#)
- [CHIR-99021](#)
- [Ciliary neurotrophic factor \(CNTF\)](#)
- [Compound E](#)
- [CX3CL1](#)
- [DAPT](#)
- [Dexamethasone](#)
- [Dibutyryl-cAMP \(db-cAMP\)](#)
- [DMEM/F12](#)
- [DMEM](#)
- [Dorsomorphin](#)
- [Dorsomorphin homologue 1 \(DMH1\)](#)
- [Embryonic stem-cell fetal bovine serum \(FBS\)](#)
- [Epidermal growth factor \(EGF\)](#)
- [Erythropoietin \(EPO\)](#)
- [Essential 8 supplement](#)
- [Essential 8 basal media](#)
- [Fetal bovine serum \(FBS\)](#)
- [Fibroblast growth factor–basic \(FGF-b\)/FGF-2](#)
- [Fibroblast growth factor 8 \(FGF-8\)](#)
- [5-Fluoro-2'-deoxyuridine](#)
- [Gentle Cell Dissociation Reagent](#)
- [Glial-derived neurotrophic factor \(GDNF\)](#)
- [GlutaMAX-I](#)
- [Heparin](#)







## Select a culture reagent to go to its screen. Culture Reagents: I–P

- [Insulin](#)
- [Insulin-like growth factor 1 \(IGF-1\)](#)
- [Insulin-transferrin-sodium selenite supplement \(ITS\)](#)
- [Interleukin-3 \(IL-3\)](#)
- [Interleukin-6 \(IL-6\)](#)
- [Interleukin-34 \(IL-34\)](#)
- [Knockout DMEM/F12](#)
- [Knockout serum replacement](#)
- [Laminin \(Invitrogen\)](#)
- [Laminin \(Sigma\)](#)
- [L-Ascorbic acid 2-phosphate \(AA2P\)](#)
- [LDN193189](#)
- [L-Glutamine](#)
- [Lithium chloride \(LiCl\)](#)
- [Macrophage colony stimulating factor \(M-CSF\)](#)
- [Matrigel growth factor reduced \(GFR\) basement membrane matrix](#)
- [Matrigel human embryonic stem cell \(hESC\)–qualified matrix](#)
- [MEM nonessential amino acid \(NEAA\) solution](#)
- [2-mercaptoethanol \( \$\beta\$ -mercaptoethanol \[BME\]; Gibco #21985023\)](#)
- [2-mercaptoethanol \( \$\beta\$ -mercaptoethanol \[BME\]; Merck #8057400005\)](#)
- [Mitomycin C](#)
- [mTeSR1 5x supplement](#)
- [mTeSR1 basal media](#)
- [N-2 supplement](#)
- [\$\beta\$ -Nerve growth factor \( \$\beta\$ -NGF\)](#)
- [Neurobasal \(NB\) media](#)
- [Neurotrophin-3 \(NT3\)](#)
- [Noggin](#)
- [Penicillin-Streptomycin](#)
- [Platelet-derived growth factor–AA \(PDGF-AA\)](#)
- [Phosphate-buffered saline \(PBS\)](#)
- [Poly-L-ornithine \(PO\)](#)
- [Polyvinyl alcohol \(PVA\)](#)
- [Purmorphamine](#)



# Materials



## Select a culture reagent to go to its screen. Culture Reagents: Q–Z

- [Retinoic acid \(RA\)](#)
- [SB431542](#)
- [Smoothened Agonist \(SAG\)](#)
- [Sodium butyrate](#)
- [Sonic hedgehog \(SHH\)](#)
- [Sonic hedgehog \(SHH; C24II\)](#)
- [Stem cell factor \(SCF\)](#)
- [Stemline II hematopoietic stem cell expansion medium](#)
- [StemPro Accutase Cell Dissociation Reagent](#)
- [Thiazovivin](#)
- [1-Thioglycerol](#)
- [Thrombopoietin \(TPO\)](#)
- [Transforming growth factor- \$\beta\$ 3 \(TGF- \$\beta\$ 3\)](#)
- [3,3',5-Triiodo-L-thyronine \(T3\)](#)
- [Valproic acid \(VPA\)](#)
- [XAV939](#)
- [Y-27632](#)





# Culture Reagents: A–H



# Activin A



## Peprotech #120-14

- Regulates cell proliferation and differentiation and promotes neuronal survival



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	26.0 kDa		<ul style="list-style-type: none"><li>Powder; briefly centrifuge vial before opening</li></ul>	-20°C
<b>Stock solution</b>	12.5 µg/mL		<ul style="list-style-type: none"><li>Dissolve 100 ug in 8 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	Microglia differentiation	12.5 ng/mL	<ul style="list-style-type: none"><li>Dilute stock solution 1:1000</li></ul>	4°C



# Antibiotic-Antimycotic



## Gibco #15240-062

- Contains the antibiotics penicillin and streptomycin and the antifungal agent amphotericin B



### MW/Concentration



### Preparation instructions



### Storage



### Important information

**As supplied/ stock solution**

100x

- Solution

-20°C  
for up to 1 year from  
date of manufacture

**Final solution (media)**

iPSC culture;  
DA neuron, cortical  
neuron, and motor  
neuron differentiation

1x

- Dilute stock solution  
1:100

4°C

**Final solution (coating)**

Laminin solution

1x

- Dilute stock solution  
1:100

4°C



# Ascorbic acid (AA)



## Sigma #A5960

- An enzymatic cofactor and antioxidant



### MW/Concentration

### Preparation instructions

### Storage

### Important information

#### As supplied

176.12 g/mol

- Powder

Room temperature

- May darken in storage

#### Stock solution

200 mM

- Dissolve 176.12 mg in 5 mL sterile ddH<sub>2</sub>O
- Filter solution
- Prepare 200- $\mu$ L aliquots in 0.7-mL tubes

-80°C

#### Final solution (media)

DA neuron or cortical neuron differentiation

200  $\mu$ M

- Dilute stock solution 1:1000

4°C

Motor neuron differentiation, midbrain organoid culture

100  $\mu$ M

- Dilute stock solution 1:2000

4°C



# B-27 supplement



## Gibco #17504044

- Supplement that supports the neuronal growth and viability



### MW/Concentration



### Preparation instructions



### Storage



### Important information

#### As supplied/ stock solution

50x

- Solution

-20°C  
for 12 months from date  
of manufacture

- Protect from light

#### Final solution (media)

DA neuron, cortical  
neuron, and  
astrocyte  
differentiation

1x

- Dilute stock solution 1:50

4°C

- Protect from light

Motor neuron  
differentiation,  
cerebral organoid  
culture

0.5x

- Dilute stock solution  
1:100

4°C

- Protect from light



# B-27 supplement, minus Vitamin A



## Gibco #12587010

- Supplement that is ideal for the cultivation of neural progenitor and stem cells, either as organoids in suspension or in adherent monolayer culture, without inducing differentiation



### MW/Concentration



### Preparation instructions



### Storage



### Important information

**As supplied/ stock solution**

50x

- Solution

-20°C  
for 12 months from date of manufacture

- Protect from light

**Final solution (media)**

Forebrain and midbrain organoid culture

1x

- Dilute stock solution 1:50

4°C

- Protect from light

Cerebral organoid culture

0.5x

- Dilute stock solution 1:100

4°C

- Protect from light





# Biotin



## Sigma #B4639

- Essential vitamin that is important for amino acid and energy metabolism and fatty acid synthesis



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	244.31 g/mol		<ul style="list-style-type: none"><li>• Powder</li></ul>	4°C
<b>Stock solution</b>	100 µg/mL		<ul style="list-style-type: none"><li>• Dissolve 10 mg in 1 mL 1N NaOH</li><li>• Add 100 µL of solution to 9.9 mL 1x PBS</li><li>• Filter solution</li><li>• Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	Oligodendrocyte differentiation	100 ng/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C



# Bone morphogenetic protein 4 (BMP4)



## GenScript #Z02750

- Regulates growth, differentiation, chemotaxis, and apoptosis of various cell types



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>		~13.3 kDa	<ul style="list-style-type: none"><li>Powder; briefly centrifuge vial before opening</li></ul>	-20°C	<ul style="list-style-type: none"><li>Stable at 4°C but should be kept at -20°C for long-term storage</li></ul>
<b>Stock solution</b>		50 µg/mL	<ul style="list-style-type: none"><li>Dissolve 100 mg in 2 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C	<ul style="list-style-type: none"><li>Stable at 4°C for up to 1 week</li><li>Avoid repeated freeze-thaw cycles</li></ul>
<b>Final solution (media)</b>	Microglia differentiation	50 ng/mL	<ul style="list-style-type: none"><li>Dilute stock solution 1:1000</li></ul>	4°C	



# Bone morphogenetic protein 9 (BMP9)/GDF-2



## Peprtech #120-07

- Regulates cell differentiation and survival



### MW/Concentration



### Preparation instructions



### Storage



### Important information

	MW/Concentration	Preparation instructions	Storage	Important information
<b>As supplied</b>	24.1 kDa	<ul style="list-style-type: none"><li>Powder; briefly centrifuge vial before opening</li></ul>	-20°C	
<b>Stock solution</b>	10 µg/mL	<ul style="list-style-type: none"><li>Dissolve 200 µg in 20 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C	
<b>Final solution (media)</b>	Cholinergic neuron differentiation 10 ng/mL	<ul style="list-style-type: none"><li>Dilute stock solution 1:1000</li></ul>	4°C	



# Bovine serum albumin (BSA)



## Invitrogen #15260-037

- Supports growth of human hematopoietic progenitor cells in serum-free media formulations



### MW/Concentration



### Preparation instructions



### Storage



### Important information

#### As supplied/ stock solution

75 mg/mL in PBS

- Liquid

-20°C

#### Final solution (media)

DA NPC and  
cortical NPC  
induction

1 mg/mL

- Dilute stock solution 1:75

4°C



# Brain-derived neurotrophic factor (BDNF)



## Peprotech #450-02

- Neurotrophic growth factor that supports neuronal proliferation and survival



### MW/Concentration



### Preparation instructions



### Storage



### Important information

**As supplied**

27.0 kDa

- Powder; briefly centrifuge vial before opening

-20°C

**Stock solution**

20 µg/mL

- Dissolve 500 µg in 25 mL sterile ddH<sub>2</sub>O containing 0.1% BSA
- Prepare 100-µL aliquots in 0.7-mL tubes

-80°C

**Final solution (media)**

DA and cortical neuron differentiation, forebrain organoid culture

20 ng/mL

- Dilute stock solution 1:1000

4°C

Motor neuron differentiation, midbrain organoid culture

10 ng/mL

- Dilute stock solution 1:2000

4°C



# CD200



## Novoprotein #C311

- Immunoregulatory protein



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	30 kDa	<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C	<ul style="list-style-type: none"><li>• Stable at room temperature for up to 3 weeks</li></ul>
<b>Stock solution</b>	100 µg/mL	<ul style="list-style-type: none"><li>• Dissolve 200 µg in 2 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>• Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C for up to 3 months	<ul style="list-style-type: none"><li>• Do not mix by vortex or pipetting</li><li>• Stable at 4°C for 2 to 7 days</li><li>• Minimize freeze-thaw cycles</li></ul>
<b>Final solution (media)</b>	Microglia differentiation 100 ng/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C	



# CHIR-99021



## Selleckchem #S2924

- Inhibits GSK-3 $\alpha/\beta$  leading to activation of GSK and insulin signaling



### MW/Concentration



### Preparation instructions



### Storage



### Important information

**As supplied**

501.4 g/mol

- Powder; briefly centrifuge vial before opening

-20°C  
for up to 3 years

**Stock solution**

3 mM

- Dissolve 25 mg in 16.607 mL DMSO
- Prepare 200- $\mu$ L aliquots in 1.5-mL tubes

-80°C  
for up to 2 years

**Final solution (media)**

DA NPC and cortical NPC induction; motor NPC induction (step 1 and 3 media)

3  $\mu$ M

- Dilute stock solution 1:1000

4°C

Motor NPC induction (step 2 media)

1  $\mu$ M

- Dilute stock solution 1:3000

4°C

Midbrain organoid culture

0.8  $\mu$ M

- Dilute stock solution 1:3750

4°C



# Ciliary neurotrophic factor (CNTF)



## Peprotech #450-13

- Neural factor that supports neuronal survival



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	22.8 kDa		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C
<b>Stock solution</b>	10 µg/mL		<ul style="list-style-type: none"><li>• Dissolve 250 µg in 25 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>• Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	Motor neuron and astrocyte differentiation	10 ng/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C





# Compound E



## STEMCELL Technologies #73954

- Inhibits  $\gamma$ -secretase and Notch processing



### MW/Concentration



### Preparation instructions



### Storage



### Important information

**As supplied**

490.5 g/mol

- Powder; briefly centrifuge vial before opening

-20°C  
for up to  
12 months

For long-term  
storage, store with a  
desiccant

- Protect from light

**Stock solution**

0.1 mM

- Dissolve 1 mg in 20.388 mL DMSO
- Prepare 100- $\mu$ L aliquots in 0.7-mL tubes

-20°C

- Avoid repeated freeze-thaw cycles
- Protect from light

**Final solution (media)**

DA neuron,  
cortical neuron,  
and motor neuron  
differentiation

0.1  $\mu$ M

- Dilute stock solution 1:1000

4°C

- Protect from light



# CX3CL1



## Peprotech #300-31

- Chemokine that chemoattracts specific cell types, including microglia cells



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	8.5 kDa	<ul style="list-style-type: none"><li>Powder; briefly centrifuge vial before opening</li></ul>	-20°C
<b>Stock solution</b>	100 µg/mL	<ul style="list-style-type: none"><li>Dissolve 500 µg in 5 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	Microglia differentiation 100 ng/mL	<ul style="list-style-type: none"><li>Dilute stock solution 1:1000</li></ul>	4°C



# DAPT



## Selleckchem #S2215

- Inhibits  $\gamma$ -secretase



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	432.46 g/mol		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C for up to 3 years
<b>Stock solution</b>	10 mM		<ul style="list-style-type: none"><li>• Dissolve 10 mg in 2.312 mL DMSO</li><li>• Prepare 100-<math>\mu</math>L aliquots in 0.7-mL tubes</li></ul>	-80°C for up to 2 years
<b>Final solution (media)</b>	DA neuron and sensory neuron differentiation	10 $\mu$ M	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C



# Dexamethasone



## Sigma #D4902

- Activates glucocorticoid receptor signaling



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	392.46 g/mol		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	4°C	
<b>Stock solution</b>	10 mM		<ul style="list-style-type: none"><li>• Dissolve 25 mg in 6.37 mL DMSO</li><li>• Prepare 20-<math>\mu</math>L aliquots in 0.7-mL tubes</li></ul>	-80°C	<ul style="list-style-type: none"><li>• Avoid repeated freeze-thaw cycles</li></ul>
<b>Final solution (media)</b>	PBMC reprogramming	1 $\mu$ M	<ul style="list-style-type: none"><li>• Dilute stock solution 1:10,000</li></ul>	4°C	



# Dibutyryl-cAMP (db-cAMP)



## Carbosynth #ND07996

- cAMP analog that activates cAMP-dependent kinases



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	491.37 g/mol	<ul style="list-style-type: none"><li>Powder; briefly centrifuge vial before opening</li></ul>	-20°C	
<b>Stock solution</b>	0.5 M	<ul style="list-style-type: none"><li>Dissolve 1 g in 4.07 mL sterile ddH<sub>2</sub>O</li><li>Prepare 200-µL aliquots in 0.7-mL tubes</li></ul>	-80°C	
<b>Final solution (media)</b>	DA neuron and cortical neuron differentiation	0.5 mM	<ul style="list-style-type: none"><li>Dilute stock solution 1:1000</li></ul>	4°C
	Midbrain organoid culture	125 µM	<ul style="list-style-type: none"><li>Dilute stock solution 1:4000</li></ul>	4°C



# Dorsomorphin



## Tocris #3093

- Inhibits AMPK and type I BMP receptors



### MW/Concentration



### Preparation instructions



### Storage



### Important information

#### As supplied

472.4 g/mol

- Powder; briefly centrifuge vial before opening

-20°C

- Product is hygroscopic. Desiccate upon arrival
- Can be stored at room temperature for up to 6 months

#### Stock solution

20 mM

- Dissolve 10 mg in 4.237 mL DMSO
- Prepare 100- $\mu$ L aliquots in 0.7-mL tubes

-80°C

- Can be stored at -20°C for up to 1 month

#### Final solution (media)

NPC induction

2  $\mu$ M

- Dilute stock solution 1:10,000

4°C

Forebrain organoid culture

10  $\mu$ M

- Dilute stock solution 1:2000

4°C



# Dorsomorphin homolog 1 (DMH1)



## Selleckchem #S7146

- Inhibits ALK2, a type 1 BMP receptor



### MW/Concentration



### Preparation instructions



### Storage



### Important instructions

#### As supplied

380.44 g/mol

- Powder; briefly centrifuge vial before opening

-20°C  
for up to 3 years

#### Stock solution

4 mM

- Dissolve 25 mg in 16.425 mL DMSO
- Prepare 100- $\mu$ L aliquots in 0.7-mL tubes

-80°C  
for up to 2 years

#### Final solution (media)

Motor neuron and astrocyte differentiation

2  $\mu$ M

- Dilute stock solution 1:2000

4°C



# DMEM/F12 media



## Gibco #10565-018

- Basal media for supporting growth of mammalian cells



### MW/Concentration



### Preparation instructions



### Storage



### Important instructions

<b>As supplied</b>		1x	<ul style="list-style-type: none"><li>Solution</li></ul>	4°C for up to 12 months from date of manufacture	<ul style="list-style-type: none"><li>Protect from light</li></ul>
<b>Final solution (coating)</b>	Laminin solution	1x	<ul style="list-style-type: none"><li>Add Antibiotic-Antimycotic and laminin as per protocol</li></ul>	4°C	<ul style="list-style-type: none"><li>Protect from light</li></ul>
<b>Final solution (media)</b>	iPSC culture; DA neuron, cortical neuron, motor neuron, and astrocyte differentiation; cerebral organoid culture	1x	<ul style="list-style-type: none"><li>Add supplements as per protocol</li></ul>	4°C	<ul style="list-style-type: none"><li>Protect from light</li></ul>





# DMEM



## Wisent Bioproducts #319-005-CL

- Basal media for supporting growth of mammalian cells



### MW/Concentration



### Preparation instructions



### Storage



### Important instructions

	MW/Concentration	Preparation instructions	Storage	Important instructions
<b>As supplied</b>	1x	<ul style="list-style-type: none"><li>Solution</li></ul>	4°C	
<b>Final solution (media)</b>	MEF culture	1x	<ul style="list-style-type: none"><li>Add supplements as per protocol</li></ul>	4°C



# Embryonic stem-cell fetal bovine serum (FBS)



## Gibco #10439024

- Sustains undifferentiated cellular morphology of embryonic stem cells



### MW/Concentration



### Preparation instructions



### Storage



### Important instructions

#### As supplied

100%

- Solution

-20°C

#### Final solution (media)

Cerebral organoid culture

3%

- Add 3 mL to 100 mL culture media

4°C



# Epidermal growth factor (EGF)



## Peprtech #AF-100-15

- Growth factor that stimulates growth of epidermal and epithelial cells



### MW/Concentration

### Preparation instructions

### Storage

### Important information

<b>As supplied</b>	6.2 kDa		<ul style="list-style-type: none"> <li>Powder; briefly centrifuge vial before opening</li> </ul>	-20°C	
<b>Stock solution</b>	10 µg/mL		<ul style="list-style-type: none"> <li>Dissolve 1 mg in 100 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li> <li>Prepare 100-µL aliquots in 0.7-mL tubes</li> </ul>	-80°C	
<b>Final solution (media)</b>	DA neuron, cortical neuron, and astrocyte differentiation*	10 ng/mL	<ul style="list-style-type: none"> <li>Dilute stock solution 1:1000</li> </ul>	4°C	
	Astrocyte differentiation*	20 ng/mL	<ul style="list-style-type: none"> <li>Dilute stock solution 1:500</li> </ul>	4°C	

\*Final solution concentration for astrocyte differentiation depends on protocol.



# Erythropoietin (EPO)



## Peprotech #100-64

- Hormone that stimulates proliferation and differentiation of erythroid progenitor cells



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	150 U/ $\mu$ g*		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C
<b>Stock solution</b>	2000 U/ $\mu$ L		<ul style="list-style-type: none"><li>• Dissolve 13.3 <math>\mu</math>g in 1 mL 1x PBS</li><li>• Prepare 50-<math>\mu</math>L aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	PBMC reprogramming	2 U/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C

\*1 U is equivalent to 6.7 ng.



# Essential 8 supplement



## Gibco #A1517001 (component #A1517101)

- Supplement for Essential 8 basal medium



### MW/Concentration



### Preparation instructions



### Storage



### Important instructions

**As supplied**

50x

- Solution

–20°C  
for up to 12 months  
from date of  
manufacture

- Store in non-frost-free freezer
- Do not re-freeze thawed solution

**Final solution (media)**

iPSC culture

1x

- Thaw Essential 8 50x supplement for about 1 hour at room temperature or overnight at 4°C. Mix gently
- Add entire bottle (10 mL) of Essential 8 50x supplement to 500 mL of Essential 8 basal media
- Prepare 40-mL aliquots in 50-mL conical tubes, seal with Parafilm, and store at –20°C

4°C  
for up to 2 weeks

–20°C  
for up to 6 months

- Do not thaw/warm supplement or complete media in a 37°C water bath
- Once supplement is thawed, use immediately or store at 4°C for up to 2 days
- Thaw complete media aliquots overnight at 4°C. Do not refreeze aliquots after thawing
- Warm complete media at room temperature







# Essential 8 basal media



## Gibco #A1517001 (component #A1516901)

- Xeno-free and feeder-free medium used for the growth and expansion of human iPSCs

	 MW/Concentration	 Preparation instructions	 Storage	 Important instructions
<b>As supplied</b>	1x	<ul style="list-style-type: none"> <li>• Solution</li> </ul>	<p>4°C for up to 12 months from date of manufacture</p>	<ul style="list-style-type: none"> <li>• Protect from light</li> </ul>
<b>Final solution (media)</b>	iPSC culture 1x	<ul style="list-style-type: none"> <li>• Thaw Essential 8 50x supplement for about 1 hour at room temperature or overnight at 4°C. Mix gently</li> <li>• Add entire bottle (10 mL) of Essential 8 50x supplement to 500 mL of Essential 8 basal media</li> <li>• Prepare 40-mL aliquots in 50-mL conical tubes, seal with Parafilm, and store at -20°C</li> </ul>	<p>4°C for up to 2 weeks</p> <p>-20°C for up to 6 months</p>	<ul style="list-style-type: none"> <li>• Do not thaw/warm supplement or complete media in a 37°C water bath</li> <li>• Once supplement is thawed, use immediately or store at 4°C for up to 2 days</li> <li>• Thaw complete media aliquots overnight at 4°C. Do not refreeze aliquots after thawing</li> <li>• Warm complete media at room temperature</li> </ul>



# Fetal bovine serum (FBS)



## ThermoFisher #12484028

- Supplement for culture media
- Cryopreservation media



### MW/Concentration



### Preparation instructions



### Storage



### Important instructions

#### As supplied/ stock solution

100%

- Solution

-20°C

#### Final solution (cryo- preservation)

iPSCs, DA NPCs,  
cortical NPCs, and  
motor NPCs

90%

- Add DMSO at 10%

Room temperature

#### Final solution (media)

MEF culture

10%

- Dilute stock solution  
1:10

4°C



# Fibroblast growth factor–basic (FGF-b)/FGF-2



## Peprtech #100-18B

- Growth factor that promotes proliferation and differentiation of a variety of cell types



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	17.2 kDa		<ul style="list-style-type: none"> <li>Powder; briefly centrifuge vial before opening</li> </ul>	−20°C
<b>Stock solution</b>	10 µg/mL		<ul style="list-style-type: none"> <li>Dissolve 1 mg in 100 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li> <li>Prepare 100-µL aliquots in 0.7-mL tubes</li> </ul>	−80°C
<b>Final solution (media)</b>	PBMC reprogramming; DA neuron, cortical neuron, and astrocyte differentiation*	10 ng/mL	<ul style="list-style-type: none"> <li>Dilute stock solution 1:1000</li> </ul>	4°C
	Astrocyte differentiation*	20 ng/mL	<ul style="list-style-type: none"> <li>Dilute stock solution 1:500</li> </ul>	4°C
	Cerebral organoid culture	4 ng/mL	<ul style="list-style-type: none"> <li>Dilute stock solution 1:2500</li> </ul>	4°C

\*Final solution concentration for astrocyte differentiation depends on protocol.





# Fibroblast growth factor 8 (FGF-8)



## Peprotech #100-25

- Growth factor that promotes cell proliferation and differentiation



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	22.5 kDa		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C
<b>Stock solution</b>	100 µg/mL		<ul style="list-style-type: none"><li>• Dissolve 500 µg in 5 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>• Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	NPC induction and midbrain organoid culture	100 ng/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C



# 5-Fluoro-2'-deoxyuridine



## Sigma #F0503

- Antineoplastic agent that inhibits thymidylate synthetase



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	246.19 g/mol		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	Room temperature	
<b>Stock solution</b>	100 mM		<ul style="list-style-type: none"><li>• Dissolve 100 mg in 406 <math>\mu</math>L sterile ddH<sub>2</sub>O</li><li>• Prepare 20-<math>\mu</math>L aliquots in 0.7-mL tubes</li></ul>	-80°C	<ul style="list-style-type: none"><li>• Stable at 4°C for up to 2 weeks</li><li>• Protect from light</li></ul>
<b>Final solution (media)</b>	DA neuron, cortical neuron, and motor neuron differentiation	1 $\mu$ M	<ul style="list-style-type: none"><li>• Dilute stock solution 1:100,000</li></ul>	4°C	<ul style="list-style-type: none"><li>• Protect from light</li></ul>



# Gentle Cell Dissociation Reagent



## STEMCELL Technologies #07174

- Enzyme-free reagent for dissociation of human embryonic stem cells or human (iPSCs) cells into cell aggregates for routine passaging or into a single-cell suspension



### MW/Concentration

**As  
supplied**

iPSCs, DA NPCs,  
cortical NPCs, and  
motor NPCs



### Preparation instructions

1x • Solution



### Storage

Room temperature



### Important instructions



# Glial-derived neurotrophic factor (GDNF)



## Peprotech #450-10

- Neurotrophic factor that promotes dopamine uptake and survival and morphological differentiation of midbrain neurons



### MW/Concentration



### Preparation instructions



### Storage



### Important information

**As supplied**

30.4 kDa

- Powder; briefly centrifuge vial before opening

-20°C

**Stock solution**

20 µg/mL

- Dissolve 500 µg in 25 mL sterile ddH<sub>2</sub>O containing 0.1% BSA
- Prepare 100-µL aliquots in 0.7-mL tubes

-80°C

**Final solution (media)**

DA neuron and cortical neuron differentiation

20 ng/mL

- Dilute stock solution 1:1000

4°C

Midbrain organoid culture

10 ng/mL

- Dilute stock solution 1:2000

4°C



# GlutaMAX-I



## Gibco #35050-061

- Supplement for adherent and suspension culture of mammalian cells that is an alternative to L-glutamine, with increased stability that improves cell health



### MW/Concentration



### Preparation instructions



### Storage



### Important information

**As supplied/ stock solution**

100x  
(200 mM in 0.85% NaCl)

- Solution

Room temperature for up to 24 months from date of manufacture

**Final solution (media)**

DA NPC and cortical NPC induction; motor neuron differentiation; midbrain and cerebral organoid culture

1x

- Dilute stock solution 1:100

4°C



# Heparin



## Sigma #H3149

- Enhances the antithrombin-mediated inactivation of proteases in the coagulation pathway



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	Mixture of polyanion chains of MW 6–30 kDa		<ul style="list-style-type: none"><li>Powder; briefly centrifuge vial before opening</li></ul>	–20°C
<b>Stock solution</b>	20 mg/mL		<ul style="list-style-type: none"><li>Dissolve 477 mg in 23.85 mL 1x PBS</li><li>Prepare 50-µL aliquots in 0.7-mL tubes</li></ul>	–80°C
<b>Final solution (media)</b>	Astrocyte differentiation	2 µg/mL	<ul style="list-style-type: none"><li>Dilute stock solution 1:10,000</li></ul>	4°C
	Midbrain and cerebral organoid culture	1 µg/mL	<ul style="list-style-type: none"><li>Dilute stock solution 1:20,000</li></ul>	4°C





# Culture Reagents: I-P



# Insulin



## Sigma #I2643

- Regulates cellular uptake, utilization, and storage of glucose, amino acids, and fatty acids and inhibits breakdown of glycogen, protein, and fat



### MW/Concentration

### Preparation instructions

### Storage

### Important information

**As supplied**

5807.57 g/mol

- Powder; briefly centrifuge vial before opening

-20°C

**Stock solution**

5 mg/mL

- Dissolve 50 mg in 10 mL sterile ddH<sub>2</sub>O containing 0.1% BSA
- Prepare 200- $\mu$ L aliquots in 0.7-mL tubes

-80°C

**Final solution (media)**

Midbrain organoid culture

2.5  $\mu$ g/mL

- Dilute stock solution 1:2000

4°C

Cerebral organoid culture

1.25  $\mu$ g/mL

- Dilute stock solution 1:4000

4°C





# Insulin-like growth factor 1 (IGF-1)



## Peprtech #100-11

- Growth factor that stimulates proliferation and differentiation of various cell types



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	7.6 kDa		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C
<b>Stock solution</b>	10 µg/mL		<ul style="list-style-type: none"><li>• Dissolve 1 mg in 100 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>• Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	PBMC reprogramming, motor neuron differentiation	10 ng/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C



# Insulin-transferrin-sodium selenite supplement (ITS)

Gibco #41400045



- Insulin promotes glucose and amino acid uptake, lipogenesis, intracellular transport, and protein and nucleic acid synthesis. Transferrin is an iron carrier and may reduce toxic levels of oxygen radicals and peroxide. Selenium, as sodium selenite, is a cofactor for glutathione peroxidase and other proteins, and an antioxidant



## MW/Concentration



## Preparation instructions



## Storage



## Important information

**As supplied/  
stock  
solution**

100x

- Solution

4°C  
for up to 18 months from  
date of manufacture

**Final  
solution  
(media)**

PBMC  
reprogramming,  
microglia  
differentiation

1x

- Dilute stock solution  
1:100

4°C

- Protect from light



# Interleukin-3 (IL-3)



## Peprtech #200-03

- Hematopoietic growth factor that promotes the survival, differentiation and proliferation of specific types of committed progenitor cells, including those of the erythroid lineage



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	15.0 kDa		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C
<b>Stock solution</b>	10 µg/mL		<ul style="list-style-type: none"><li>• Dissolve 1 mg in 100 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>• Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	PBMC reprogramming	10 ng/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C



# Interleukin-6 (IL-6)



## Peprotech #200-06

- Cytokine with diverse biological functions



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	20.9 kDa		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C
<b>Stock solution</b>	50 µg/mL		<ul style="list-style-type: none"><li>• Dissolve 500 µg in 10 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>• Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	Microglia differentiation	50 ng/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C



# Interleukin-34 (IL-34)



## Peprtech #200-34

- Ligand for colony-stimulating factor-1 receptor (CSF1R)



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>		52.5 kDa	<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C	
<b>Stock solution</b>		100 µg/mL	<ul style="list-style-type: none"><li>• Dissolve 500 µg in 5 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>• Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C	
<b>Final solution (media)</b>	Microglia differentiation	100 ng/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C	



# KnockOut serum replacement



## Gibco #10828028

- Supplement for KnockOut DMEM/F12 basal media



### MW/Concentration



### Preparation instructions



### Storage



### Important information

**As supplied**

100%

- Solution
- Prepare 10-mL aliquots in 15-mL tubes

-20°C  
for up to 18 months  
from date of  
manufacture

- Avoid repeated freeze-thaw cycles
- Can be stored at 4°C for up to 4 weeks
- Protect from light

**Final solution (media)**

PBMC reprogramming, iPSC culture, cerebral organoid culture

20%

- Thaw aliquot of KnockOut serum replacement overnight at 4°C
- Add 20 mL KnockOut serum replacement to 100 mL KnockOut DMEM/F12
- Add supplements as per protocol

4°C  
for up to 10 days

- Do not thaw/warm serum replacement or complete media in a 37°C water bath
- Warm only the volume of complete media required for that day's use
- Protect complete media from light



# KnockOut DMEM/F12



## Gibco #12660012

- Low osmolality medium without L-glutamine or HEPES buffer for the growth of human embryonic stem cells and iPSCs



### MW/Concentration



### Preparation instructions



### Storage



### Important information

		MW/Concentration	Preparation instructions	Storage	Important information
<b>As supplied</b>		1x	<ul style="list-style-type: none"> <li>• Solution</li> </ul>	4°C	
<b>Final solution (media)</b>	PBMC reprogramming, iPSC culture	1x	<ul style="list-style-type: none"> <li>• Thaw aliquot of KnockOut serum replacement overnight at 4°C</li> <li>• Add 20 mL KnockOut serum replacement to 100 mL KnockOut DMEM/F12</li> <li>• Add supplements as per protocol</li> </ul>	4°C for up to 10 days	<ul style="list-style-type: none"> <li>• Do not thaw/warm serum replacement or complete media in a 37°C water bath</li> <li>• Warm only the volume of complete media required for that day's use</li> <li>• Protect complete media from light</li> </ul>



# Laminin



## Invitrogen #23017-015

- Extracellular matrix protein that supports adhesion, proliferation, and differentiation of many cell types



### MW/Concentration



### Preparation instructions



### Storage



### Important information

#### As supplied/stock solution

0.5–2.0 mg/mL\* in 50 mM Tris-HCl (pH 7.4), 0.15 M NaCl

- Solution
- Thaw solution slowly at 4°C
- Prepare 100-µL aliquots in 0.7-mL polypropylene microcentrifuge tubes

–80°C

- If frozen solution is warmed too quickly, laminin may form a gel and cannot be reactivated for use
- Do not freeze-thaw repeatedly

#### Final solution (coating)

DA neuron, cortical neuron, motor neuron, and astrocyte differentiation

5 µg/mL

- Dilute stock solution to 5 µg/mL in DMEM/F12

4°C for up to 2 weeks

- Do not store for more than 2 weeks

#### Final solution (media)

DA neuron, cortical neuron, and motor neuron differentiation

1 µg/mL

- Dilute stock solution to 1 µg/mL

4°C for up to 2 weeks

- Do not store for more than 2 weeks

\*Verify exact concentration on tube label.





# Laminin



## Sigma #L2020

- Extracellular matrix protein that supports adhesion, proliferation, and differentiation of many cell types



### MW/Concentration



### Preparation instructions



### Storage



### Important information

#### As supplied/ stock solution

~1 mg/mL\* in  
50 mM Tris-HCl (pH 7.5),  
150 mM NaCl

- Solution
- Thaw solution slowly at 4°C
- Prepare 100-µL aliquots in 0.7-mL polypropylene microcentrifuge tubes

-80°C

- If the frozen solution is warmed too quickly, laminin may form a gel and cannot be reactivated for use

#### Final solution (coating)

DA NPC and cortical NPC induction; DA neuron, cortical neuron, motor neuron, and astrocyte differentiation

5  
µg/mL

- Dilute stock solution to 5 µg/mL in DMEM/F12

4°C  
for up to 2  
weeks

- Do not use if discoloration or web formations appear on coated surfaces

#### Final solution (media)

DA neuron, cortical neuron, and motor neuron differentiation

1  
µg/mL

- Dilute stock solution to 1 µg/mL

4°C  
for up to 2  
weeks

- Do not store for more than 2 weeks

Midbrain organoid culture

200  
ng/mL

- Dilute stock solution to 200 ng/mL

4°C  
for up to 2  
weeks

- Do not store for more than 2 weeks

\*Verify exact concentration on tube label.



# L-Ascorbic acid 2-phosphate (AA2P; sesquimagnesium salt hydrate)

Sigma #A8960



- Long-acting ascorbic acid derivative that stimulates collagen expression and formation



## MW/Concentration



## Preparation instructions



## Storage



## Important information

<b>As supplied</b>	289.54 g/mol		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C
<b>Stock solution</b>	32 mg/mL		<ul style="list-style-type: none"><li>• Dissolve 32 mg in 1 mL sterile ddH<sub>2</sub>O</li><li>• Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	4°C
<b>Final solution (media)</b>	PBMC reprogramming	64 µg/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:500</li></ul>	4°C



# LDN193189



## Sigma #SML0559

- Dorsomorphin derivative that inhibits ALK2 and ALK3



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	406.48 g/mol		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C
<b>Stock solution</b>	1 mM		<ul style="list-style-type: none"><li>• Dissolve 5 mg in 12.3 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>• Prepare 100-μL aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	Astrocyte differentiation	0.1 μM	<ul style="list-style-type: none"><li>• Dilute stock solution 1:10,000</li></ul>	4°C



# L-Glutamine



## Wisent Bioproducts #609-065-EL

- Amino acid required for cell culture involved in the formation of purine and pyrimidine nucleotides, amino sugars, glutathione, L-glutamate, other amino acids; protein synthesis; and glucose production



MW/Concentration



Preparation instructions



Storage



Important information

**As supplied/ stock solution**

200 mM

- Solution

-20°C

**Final solution (media)**

PBMC reprogramming, MEF culture

2 mM

- Dilute stock solution 1:100

4°C



# Lithium chloride (LiCl)



## Sigma #L7026

- Greatly enhances the generation of iPSCs



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	8 M in H <sub>2</sub> O		• Solution	4°C
<b>Stock solution</b>	2 M		• Dilute in sterile ddH <sub>2</sub> O to achieve 2 M • Filter solution • Prepare 100-µL aliquots in 0.7-mL tubes	4°C
<b>Final solution (media)</b>	PBMC reprogramming	2 mM	• Dilute stock solution 1:1000	4°C



# Macrophage colony stimulating factor (M-CSF)



## Peprtech #300-25

- Hematopoietic factor that regulates cellular proliferation, differentiation, and survival in blood monocytes, tissue macrophages, and their respective progenitor cells



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	36.8 kDa		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C
<b>Stock solution</b>	25 µg/mL		<ul style="list-style-type: none"><li>• Dissolve 500 µg in 20 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>• Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	Microglia differentiation	25 ng/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C



# Matrigel growth factor reduced (GFR) basement membrane matrix



## Corning #356230

- Solubilized basement membrane extract containing laminin (a major component), collagen IV, heparin sulfate proteoglycans, entactin/nidogen, and a number of growth factors



### MW/Concentration



### Preparation instructions



### Storage



### Important information

**As supplied/stock solution**

1x in DMEM,  
50 µg/mL gentamycin

- Solution

-20°C

- Do not store in frost-free freezer

**Final solution (coating)**

Midbrain and cerebral organoid culture

1x

- Thaw stock solution overnight at 4°C on ice

4°C

- Chill pipette tips at -20°C for 3 hours
- Use entire bottle after thawing



# Matrigel human embryonic stem cell (hESC)–qualified matrix



## Corning Millipore #354277

- Soluble basement membrane extract containing laminin, collagen IV, entactin and heparan sulfate, and proteoglycan



### MW/Concentration



### Preparation instructions



### Storage



### Important information

#### As supplied/ stock solution

100x\* in DMEM,  
50 µg/mL gentamycin

- Solution
- Thaw overnight at 4°C on ice
- Prepare 150-µL aliquots in pre-chilled 0.7-mL tubes on ice

- As supplied: –20°C
- Aliquots: –80°C

- Chill tubes and pipette tips at –20°C for 3 hours
- Prepare one time–use aliquots to minimize freeze-thaw cycles
- Do not store in frost-free freezer

#### Final solution (coating)

iPSC culture,  
motor neuron and  
astrocyte  
differentiation

1x

- Thaw aliquot of stock solution on ice
- Dilute stock solution 1:100

4°C

- Use immediately after preparation or store at –20°C for up to 1 week. Thaw on ice and use immediately. Do not re-freeze

\*Exact concentration on tube label.





# 2-Mercaptoethanol ( $\beta$ -mercaptoethanol [BME])



Gibco #21985023

- Reducing agent used in culture media to prevent toxic levels of oxygen radicals



## MW/Concentration

## Preparation instructions

## Storage

## Important information

### As supplied/ stock solution

55 mM in Dulbecco's phosphate-buffered saline (DPBS)

- Solution

4°C  
for 36 months from date of manufacture

### Final solution (media)

Forebrain organoid culture

0.1 nM

- Prepare serial dilutions of stock solution in sterile ddH<sub>2</sub>O (e.g. 1:1000 to 1:50,000)
- Dilute appropriate serial dilution to achieve 0.1 nM

4°C

- The 1:1000 serial dilution of the stock solution (55  $\mu$ M) may be stored at 4°C for up to 3 months when used for the same batch of organoids. Prepare fresh serial dilutions for a new batch of organoids.



# 2-Mercaptoethanol ( $\beta$ -mercaptoethanol [BME])



Merck #8057400005

- Reducing agent used in culture media to prevent toxic levels of oxygen radicals



MW/Concentration

Preparation instructions

Storage

Important information

As  
supplied/  
stock  
solution

78.12 g/mol (14.25 mM)

- Solution

Room  
temperature

Final  
solution  
(media)

Midbrain organoid  
culture

0.00035%

- Dilute stock solution to  
0.00035%

4°C

Cerebral organoid  
culture

0.00035%  
or  
0.0007%

- Dilute stock solution to  
0.00035% or 0.0007%, as  
per protocol

4°C



# MEM nonessential amino acid (NEAA) solution



## Wisent #321-011-EL

- Supplement to increase cell growth and viability



MW/Concentration



Preparation instructions



Storage



Important information

**As supplied/  
stock  
solution**

100x

- Solution

4°C

**Final  
solution  
(media)**

PBMC reprogramming;  
DA NPC and cortical NPC  
induction; DA neuron,  
cortical neuron, and  
astrocyte differentiation;  
midbrain and cerebral  
organoid culture\*

1x

- Dilute stock solution  
1:100

4°C

Cerebral organoid  
culture\*

0.5x

- Dilute stock solution  
1:200

4°C

\*Final solution concentration for cerebral organoid culture depends on cerebral organoid stage, as per protocol.



# Mitomycin C



## Sigma #M4287

- Antibiotic and double-stranded DNA alkylating agent that inhibits cell proliferation



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	334.33 g/mol	<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	4°C		
<b>Stock solution</b>	1 mg/mL	<ul style="list-style-type: none"><li>• Dissolve 2 mg in 2 mL sterile ddH<sub>2</sub>O</li><li>• Prepare 50-µL aliquots in 0.7-mL tubes</li></ul>	-80°C	<ul style="list-style-type: none"><li>• Can be stored at 4°C for up to 1 week (pH 6–9)</li><li>• If a precipitate forms, prepare a fresh stock solution</li><li>• Protect from light</li></ul>	
<b>Final solution (media)</b>	DA neuron and cortical neuron differentiation	1 µg/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C	<ul style="list-style-type: none"><li>• Protect from light</li></ul>







# mTeSR1 5x supplement



## STEMCELL Technologies #85850 (component #85852)

- Supplement for mTeSR1 basal media

	 MW/Concentration	 Preparation instructions	 Storage	 Important instructions
<b>As supplied</b>	5x	<ul style="list-style-type: none"> <li>• Solution</li> </ul>	<p>–20°C for up to 3 months</p>	
<b>Final solution (media)</b>	iPSC culture 1x	<ul style="list-style-type: none"> <li>• Thaw mTeSR1 5x supplement overnight at 4°C. Mix thoroughly</li> <li>• Add 100 mL mTeSR1 5x supplement to 400 mL mTeSR1 basal media to obtain 500 mL of complete media. Mix well</li> <li>• Prepare 40-mL aliquots in 50-mL conical tubes, seal with Parafilm, and store at –20°C</li> </ul>	<p>4°C for up to 2 weeks</p> <p>–20°C for up to 6 months</p>	<ul style="list-style-type: none"> <li>• Do not thaw/warm supplement or complete media in a 37°C water bath</li> <li>• Once supplement is thawed, use immediately or store at 4°C for up to 2 days</li> <li>• Thaw complete media aliquots overnight at 4°C. Do not refreeze aliquots after thawing</li> <li>• Warm complete media at room temperature</li> </ul>







# mTeSR1 basal media



## STEMCELL Technologies #85850 (component #85851)

- Feeder-free cell culture media for human embryonic stem cells and iPSCs

	 MW/Concentration	 Preparation instructions	 Storage	 Important instructions
<b>As supplied</b>	1x	<ul style="list-style-type: none"> <li>Solution</li> </ul>	4°C	
<b>Final solution (media)</b>	iPSC culture 1x	<ul style="list-style-type: none"> <li>Thaw mTeSR1 5x supplement overnight at 4°C. Mix thoroughly</li> <li>Add 100 mL mTeSR1 5x supplement to 400 mL mTeSR1 basal media to obtain 500 mL of complete media. Mix well</li> <li>Prepare 40-mL aliquots in 50-mL conical tubes, seal with Parafilm, and store at -20°C</li> </ul>	4°C for up to 2 weeks  -20°C for up to 6 months	<ul style="list-style-type: none"> <li>Do not thaw/warm supplement or complete media in a 37°C water bath</li> <li>Once supplement is thawed, use immediately or store at 4°C for up to 2 days</li> <li>Thaw complete media aliquots overnight at 4°C. Do not refreeze aliquots after thawing</li> <li>Warm complete media at room temperature</li> </ul>



# N-2 supplement



## Gibco #17502048

- Supplement that promotes neuronal cell growth



### MW/Concentration



### Preparation instructions



### Storage



### Important information

**As supplied/  
stock  
solution**

100x

- Solution

-20°C

- Protect from light

**Final  
solution  
(media)**

DA neuron, cortical neuron, and astrocyte differentiation; midbrain organoid culture

1x

- Dilute stock solution 1:100

4°C

- Protect from light

Motor neuron differentiation, cerebral organoid culture

0.5x

- Dilute stock solution 1:200

4°C

- Protect from light



# $\beta$ -Nerve growth factor ( $\beta$ -NGF)



## Peprotech #450-01

- Neurotrophic factor that is crucial for development and preservation of sensory and sympathetic nervous systems



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	13.5 kDa		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C
<b>Stock solution</b>	100 $\mu$ g/mL		<ul style="list-style-type: none"><li>• Dissolve 500 <math>\mu</math>g in 5 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>• Prepare 100-<math>\mu</math>L aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	Cholinergic neuron differentiation	100 ng/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C





# Neurobasal (NB) media



## Life Technologies #21103-049

- Basal media for long-term maintenance and maturation of pure pre-natal and embryonic neuronal cell populations



### MW/Concentration



### Preparation instructions



### Storage



### Important information

**As supplied**

1x

- Solution

4°C

**Final solution (media)**

DA neuron, cortical neuron, and motor neuron differentiation; cerebral organoid culture

1x

- Add supplements as per protocol

4°C



# Neurotrophin-3 (NT3)



## Peprtech #450-03

- Neurotrophic factor that promotes the growth and survival of nerve and glial cells



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	13.6 kDa		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C
<b>Stock solution</b>	10 µg/mL		<ul style="list-style-type: none"><li>• Dissolve 50 µg in 5 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>• Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	Sensory neuron and oligodendrocyte differentiation	10 ng/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C
	Forebrain organoid culture	20 ng/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:500</li></ul>	4°C



# Noggin



## Peprotech #120-10C

- Inhibits TGF- $\beta$  ligands during developmental processes



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	46 kDa		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C
<b>Stock solution</b>	200 $\mu$ g/mL		<ul style="list-style-type: none"><li>• Dissolve 500 <math>\mu</math>g in 2.5 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>• Prepare 100-<math>\mu</math>L aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	NPC induction, midbrain organoid culture	200 ng/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C



# Penicillin-Streptomycin



## Wisent #450-200-EL

- Contains the antibiotics penicillin and streptomycin



### MW/Concentration



### Preparation instructions



### Storage



### Important information

**As  
supplied/  
stock  
solution**

100X

- Solution

-20°C

**Final  
solution  
(media)**

PBMC reprogramming,  
MEF culture, midbrain and  
cerebral organoid culture

1x

- Dilute stock solution  
1:100

4°C



# Platelet-derived growth factor-AA (PDGF-AA)



## Peprotech #100-13A

- Mitogen involved in a number of biological processes, including embryonic neuron development



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	28.5 kDa		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C
<b>Stock solution</b>	10 µg/mL		<ul style="list-style-type: none"><li>• Dissolve 50 µg in 5 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>• Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	Oligodendrocyte differentiation	10 ng/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C



# Phosphate-buffered saline (PBS)



## Wisent #311-010-CL

- Balanced salt solution used for a variety of cell culture applications, such as washing cells and preparing reagents



MW/Concentration

As supplied

1x



Preparation instructions

- Solution



Storage

Room temperature



Important information



# 2-Phospho-L-ascorbic acid (AA2P) trisodium salt



Sigma #49752

- Long-acting ascorbic acid derivative that stimulates collagen expression and formation



MW/Concentration



Preparation instructions



Storage



Important information

<b>As supplied</b>	322.05 g/mol		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C
<b>Stock solution</b>	50 mg/mL		<ul style="list-style-type: none"><li>• Dissolve 500 mg in 10 mL sterile ddH<sub>2</sub>O</li><li>• Filter solution</li><li>• Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	PBMC reprogramming	50 µg/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C



# Poly-L-ornithine (PO)



## Sigma #P3655

- Synthetic amino acid polymer coating that acts as a charge modifier



### MW/Concentration

### Preparation instructions

### Storage

### Important information

<b>As supplied</b>	30–70 kDa	<ul style="list-style-type: none"><li>• Powder</li></ul>	–20°C	
<b>Stock solution</b>	1 mg/mL	<ul style="list-style-type: none"><li>• Dissolve 100 mg in 100 mL sterile 1x PBS</li><li>• Prepare 10-mL aliquots in 15-mL tubes</li></ul>	–20°C	
<b>Final solution (coating)</b>	DA neuron, cortical neuron, motor neuron, and astrocyte differentiation 10 µg/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:100 in 1x PBS</li></ul>	4°C for up to 2 months	<ul style="list-style-type: none"><li>• Ensure transfer of all stock solution by rinsing tube with 1x PBS twice</li></ul>





# Polyvinyl alcohol (PVA)



## Sigma #P8136

- Hydrophilic linear polymer which forms copolymers of vinyl alcohol and vinyl acetate



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	28.5 kDa	<ul style="list-style-type: none"><li>Powder</li></ul>	-20°C
<b>Stock solution</b>	10 mg/mL	<ul style="list-style-type: none"><li>Dissolve 100 mg in 10 mL ddH<sub>2</sub>O</li><li>Filter solution</li><li>Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (coating)</b>	Microglia differentiation 10 µg/mL	<ul style="list-style-type: none"><li>Dilute stock solution 1:1000</li></ul>	4°C



# Purmorphamine



## Sigma #SML-0868

- Binds to the Smo receptor and activates Hedgehog signaling



### MW/Concentration



### Preparation instructions



### Storage



### Important information

	MW/Concentration	Preparation instructions	Storage	Important information
<b>As supplied</b>	520.6 g/mol	<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C	
<b>Stock solution</b>	2 mM	<ul style="list-style-type: none"><li>• Dissolve 5 mg in 4.802 mL DMSO</li><li>• Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C	<ul style="list-style-type: none"><li>• Prepare as accurately as possible because the working concentration range is very narrow</li></ul>
<b>Final solution (media)</b>	DA neuron differentiation	2 µM	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C
	Motor NPC induction (step 2 and 3 media)	0.5 µM	<ul style="list-style-type: none"><li>• Dilute stock solution 1:4000</li></ul>	4°C
	Motor neuron differentiation	0.1 µM	<ul style="list-style-type: none"><li>• Dilute stock solution 1:20,000</li></ul>	4°C





# Culture Reagents: Q–Z



# Retinoic acid (RA)



## Sigma #R2625

- Activates transcription factors that regulate cell growth and differentiation



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>		300.44 g/mol	<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C for up to 1 year	<ul style="list-style-type: none"><li>• Stable if unopened in ampule</li><li>• Protect from light and air</li></ul>
<b>Stock solution</b>		1 mM	<ul style="list-style-type: none"><li>• Dissolve 300 mg in 10 mL DMSO</li><li>• Dilute solution 1:10 in EtOH</li><li>• Prepare 20-µL aliquots in 0.7-mL tubes</li></ul>	-80°C	<ul style="list-style-type: none"><li>• Use all of the powder immediately after opening the ampule</li><li>• Protect from light and air</li></ul>
<b>Final solution (media)</b>	Motor neuron differentiation	0.5 µM	<ul style="list-style-type: none"><li>• Dilute stock solution 1:2000</li></ul>	4°C	<ul style="list-style-type: none"><li>• Always use a new aliquot of stock solution to prepare the final solution</li><li>• Protect from light and air</li></ul>
	Motor NPC induction (step 2 and 3 media)	0.1 µM	<ul style="list-style-type: none"><li>• Dilute stock solution 1:10,000</li></ul>	4°C	



# SB431542



## Selleckchem #S1067

- Inhibits ALK5 and the TGF- $\beta$ /Activin/NODAL pathway



### MW/Concentration

### Preparation instructions

### Storage

### Important information

**As supplied**

384.4 g/mol

- Powder; briefly centrifuge vial before opening

-20°C  
for up to  
3 years

**Stock solution**

10 mM

- Dissolve 10 mg in 2.604 mL DMSO
- Prepare 100- $\mu$ L aliquots in 0.7-mL tubes

-80°C  
for up to  
2 years

**Final solution (media)**

DA NPC and cortical NPC induction, astrocyte differentiation\*, midbrain organoid culture

10  $\mu$ M

- Dilute stock solution 1:1000

4°C

Motor neuron and astrocyte differentiation\*

2  $\mu$ M

- Dilute stock solution 1:5000

4°C

\*Final solution concentration for astrocyte differentiation depends on protocol.



# Smoothened Agonist (SAG)



## Millipore #566660

- Activates Smo signalling, Smo internalization, and Hedgehog signaling



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>		599 g/mol	<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C	<ul style="list-style-type: none"><li>• Can be stored at 4°C</li><li>• Protect from light</li></ul>
<b>Stock solution</b>		1 mM	<ul style="list-style-type: none"><li>• Dissolve 1 mg in 1.669 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>• Prepare 100-μL aliquots in 0.7-mL tubes</li></ul>	-80°C	<ul style="list-style-type: none"><li>• Can be stored at -20°C for up to 6 months</li><li>• Protect from light</li></ul>
<b>Final solution (media)</b>	NPC induction	0.25 μM	<ul style="list-style-type: none"><li>• Dilute stock solution 1:4000</li></ul>	4°C	<ul style="list-style-type: none"><li>• Protect from light</li></ul>



# Sodium butyrate



## Sigma #B5887

- Inhibits histone deacetylases



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	110 g/mol		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	Room temperature
<b>Stock solution</b>	125 mM		<ul style="list-style-type: none"><li>• Dissolve 250 mg in 18.18 mL sterile ddH<sub>2</sub>O</li><li>• Prepare 100-<math>\mu</math>L aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	PBMC reprogramming	250 $\mu$ M	<ul style="list-style-type: none"><li>• Dilute stock solution 1:500</li></ul>	4°C



# Sonic hedgehog (SHH)



## Peprtech #100-45

- Morphogen that activates Hedgehog signaling and is involved in central nervous system patterning during development



MW/Concentration



Preparation instructions



Storage



Important information

**As supplied**

20.0 kDa

- Powder; briefly centrifuge vial before opening

-20°C

**Stock solution**

200 µg/mL

- Dissolve 500 µg in 2.5 mL sterile ddH<sub>2</sub>O containing 0.1% BSA
- Prepare 100-µL aliquots in 0.7-mL tubes

-80°C

**Final solution (media)**

NPC induction

200 ng/mL

- Dilute stock solution 1:1000

4°C

Midbrain organoid culture

100 ng/mL

- Dilute stock solution 1:2000

4°C





# Sonic hedgehog (SHH; C24II)



## GenScript #Z03067

- Morphogen that activates Hedgehog signaling and is involved in central nervous system patterning during development



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	20.0 kDa	<ul style="list-style-type: none"><li>Powder; briefly centrifuge vial before opening</li></ul>	-80°C for up to 6 months from date of receipt	
<b>Stock solution</b>	200 µg/mL	<ul style="list-style-type: none"><li>Dissolve 500 µg in 2.5 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C	<ul style="list-style-type: none"><li>Can be stored at 4°C for up to 2 weeks or at -20°C for up to 3 months</li></ul>
<b>Final solution (media)</b>	NPC induction 200 ng/mL	<ul style="list-style-type: none"><li>Dilute stock solution 1:1000</li></ul>	4°C	



# Stem cell factor (SCF)



## Peprotech #300-07

- Hematopoietic growth factor that is essential for the survival, proliferation and differentiation of hematopoietic cells committed to the melanocyte and germ cell lineages



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	18.4 kDa		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C
<b>Stock solution</b>	10 µg/mL		<ul style="list-style-type: none"><li>• Dissolve 50 µg in 5 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>• Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	PBMC reprogramming	100 ng/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:100</li></ul>	4°C



# Stemline II hematopoietic stem cell expansion medium

Sigma #S0192



- Hematopoietic stem cell expansion media for differentiated and undifferentiated cells



MW/Concentration



Preparation instructions



Storage



Important information

**As supplied**

1x

- Solution

4°C

- Protect from light

**Final solution  
(media)**

1x

- Add supplements as per protocol

4°C

- Protect from light



# StemPro Accutase Cell Dissociation Reagent



## ThermoFisher #A1110501

- Cell detachment solution of proteolytic and collagenolytic enzymes



### MW/Concentration

**As supplied**

DA NPCs and cortical NPCs, midbrain organoids

1x



### Preparation instructions

- Solution



### Storage

- Upon arrival:  $-20^{\circ}\text{C}$
- After thawing:  $4^{\circ}\text{C}$  for up to 2 years



### Important information

- Multiple freeze-thaw cycles are not recommended
- Do not store at room temperature
- Protect from light



# Thiazovivin



## Tocris #3845

- Inhibits ROCK and the RHO/ROCK pathway



### MW/Concentration



### Preparation instructions



### Storage



### Important information

**As supplied**

311 g/mol

- Powder; briefly centrifuge vial before opening

-20°C  
for up to  
6 months

**Stock solution**

2 mM

- Dissolve 10 mg in 16.367 mL DMSO
- Prepare 100-µL aliquots in 0.7-mL tubes

-20°C  
for up to  
1 month

**Final solution (media)**

iPSC culture; DA neuron, cortical neuron, and motor neuron differentiation

2 µM

- Dilute stock solution 1:1000

4°C



# 1-Thioglycerol



## Sigma #M1753

- Stimulates cell proliferation



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	1.25 g/mL (11.57 M)		<ul style="list-style-type: none"><li>Solution</li></ul>	4°C
<b>Stock solution</b>	0.4 M		<ul style="list-style-type: none"><li>Dilute in sterile ddH<sub>2</sub>O to achieve 0.4 M</li><li>Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	4°C
<b>Final solution (media)</b>	PBMC reprogramming	0.4 mM	<ul style="list-style-type: none"><li>Dilute stock solution 1:1000</li></ul>	4°C



# Thrombopoietin (TPO)



## Peprtech #300-18

- Growth factor that stimulates the proliferation and maturation of megakaryocytes, and promotes increased circulating levels of platelets



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	18.4 kDa		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C
<b>Stock solution</b>	50 µg/mL		<ul style="list-style-type: none"><li>• Dissolve 250 µg in 5 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>• Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	Microglia differentiation	50 ng/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C



# Transforming growth factor $\beta 3$ (TGF- $\beta 3$ )



## Peprtech #100-36E

- Cytokine that regulates cell proliferation, growth, differentiation and motility, as well as synthesis and deposition of the extracellular matrix



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	25.0 kDa		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C
<b>Stock solution</b>	1 $\mu$ g/mL		<ul style="list-style-type: none"><li>• Dissolve 50 <math>\mu</math>g in 50 mL sterile ddH<sub>2</sub>O containing 0.1% BSA</li><li>• Prepare 100-<math>\mu</math>L aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	DA neuron and cortical neuron differentiation	1 ng/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C





# 3,3',5-Triiodo-L-thyronine (T3)



## Sigma #T6397

- Thyroid hormone that regulates cell differentiation and protein expression



### MW/Concentration

### Preparation instructions

### Storage

### Important information

<b>As supplied</b>	672.96 g/mol		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C	
<b>Stock solution</b>	60 µg/mL		<ul style="list-style-type: none"><li>• Dissolve 6 mg in 10 mL 0.1N NaOH</li><li>• Add 1 mL of solution to 9 mL 1x PBS</li><li>• Filter solution</li><li>• Prepare 100-µL aliquots in 0.7-mL tubes</li></ul>	-80°C	<ul style="list-style-type: none"><li>• Can be stored at 4°C for up to 1 month or at -20°C</li><li>• Avoid repeated freezing and thawing of aliquots</li></ul>
<b>Final solution (media)</b>	Oligodendrocyte differentiation	60 ng/mL	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C	



# Valproic acid (VPA)



## Sigma #P4543

- Anticonvulsant; inhibits histone deacetylases and enzymes involved in GABA metabolism



### MW/Concentration

### Preparation instructions

### Storage

### Important information

<b>As supplied</b>	166.2 g/mol	<ul style="list-style-type: none"><li>• Powder</li></ul>	Room temperature
<b>Stock solution</b>	0.5 M	<ul style="list-style-type: none"><li>• Dissolve 831 mg in 10 mL ddH<sub>2</sub>O</li><li>• Filter solution</li><li>• Prepare 100-<math>\mu</math>L aliquots in 0.7-mL tubes</li></ul>	-80°C
<b>Final solution (media)</b>	DA neuron, cortical neuron, and motor neuron differentiation	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C
	Astrocyte differentiation	<ul style="list-style-type: none"><li>• Dilute stock solution 1:50,000</li></ul>	4°C



# XAV939



## Sigma #X3004

- Inhibits Tankyrase and Wnt/ $\beta$ -catenin signaling



### MW/Concentration



### Preparation instructions



### Storage



### Important information

<b>As supplied</b>	312 g/mol		<ul style="list-style-type: none"><li>• Powder; briefly centrifuge vial before opening</li></ul>	-20°C	<ul style="list-style-type: none"><li>• Can be stored at room temperature</li></ul>
<b>Stock solution</b>	10 mM		<ul style="list-style-type: none"><li>• Dissolve 25 mg in 8.012 mL DMSO</li><li>• Prepare 100-<math>\mu</math>L aliquots in 0.7-mL tubes</li></ul>	-80°C	
<b>Final solution (media)</b>	GABAergic neuron differentiation	1 $\mu$ M	<ul style="list-style-type: none"><li>• Dilute stock solution 1:1000</li></ul>	4°C	



# Y-27632



## Selleckchem #S1049

- Inhibits ROCK



### MW/Concentration



### Preparation instructions



### Storage



### Important information

**As supplied**

320.26 g/mol

- Powder; briefly centrifuge vial before opening

–20°C  
for up to  
3 years

**Stock solution**

10 mM

- Dissolve 50 mg in 15.612 mL sterile ddH<sub>2</sub>O
- Prepare 200-µL aliquots in 0.7-mL tubes

–80°C  
for up to  
2 years

**Final solution (media)**

iPSC culture; DA neuron, cortical neuron, and motor neuron differentiation; midbrain organoid culture

10 µM

- Dilute stock solution 1:1000

4°C

Cerebral organoid culture

50 µM

- Dilute stock solution 1:200

4°C





**You have reached the end of the protocol.**

