

Username : Quentin Laurent

Confidential

**QL13 QL13-11**

Done

**Title:** QL13-11**Goal:** PTO Cy5 uptake, 500 nM, activation + activation followed by dialysis with 5 compounds**Started on:** March 18, 2021**Modified on:** May 12, 2021*Closed date***Project:** Phosphorothioates**Type:** PTO Cy5 Inhibition

HeLa Kyoto cells were seeded at  $8 \times 10^4$  cells/mL in FluoroBrite DMEM + 10% FBS on  $\mu$ -Plates 96-well ibiTreat sterile and kept at 37 °C with 5% CO<sub>2</sub> overnight. The following day, cells were washed with FluoroBrite DMEM (3 × 3 mL/well) using a plate washer, keeping a final volume of 135  $\mu$ L/well.

Dilutions of Cy5-labelled phosphorothioate and inhibitor in PBS were prepared in a 96-well V-bottom plate and added to the  $\mu$ -Plate containing the cells at the desired time (15  $\mu$ L/well, 10x final concentration in PBS) to reach a final volume of 150  $\mu$ L/well, and cells were incubated for 1 h at 37 °C with 5% CO<sub>2</sub> in DMEM. Cells were washed with FluoroBrite DMEM (3 × 3 mL/well) using a plate washer, keeping a final volume of 135  $\mu$ L/well, and the corresponding solution of phosphorothioate was added (15  $\mu$ L/well, 10x final concentration in PBS) to reach a final volume of 150  $\mu$ L/well, and cells were incubated for 2 h at 37 °C with 5% CO<sub>2</sub> in DMEM.

Then, Hoechst 33342 (15  $\mu$ L/well, 100  $\mu$ g/mL in PBS) was added. After 15 min of incubation at 37 °C with 5% CO<sub>2</sub>, the medium was removed and cells were fixed with a PFA solution (3% in PBS, 80  $\mu$ L/well) for 20 min. The plate was washed with PBS (3 x 3 mL/well) and imaged using a high-content automated microscope.

## Attached files

QL13-01-Plates.key

QL13-01-Activation.pzfx

QL13-01-Inhibition.pzfx