



Continued from page:

96-well plate inhibition


1 h pre-incubation with inhibitors, 2 h incubation PTO

Cells seeded at  $8 \cdot 10^4$ , incubated 24 h at 37 °C

Stock sol<sup>o</sup> 10 mM O=C(O)c1ccccc1 → ETP-S<sub>4</sub> → 100 μM

Stock sol<sup>o</sup> 10 mM O=C(=O)c1ccccc1 → 1 mM

Cells + 135 μL medium, add 15 μL 10x inhibitor



1 h incubation inhibitors, wash with plate washer

Add 15 μL 5 μM PS/PO in each well

Stock sol<sup>o</sup> 100 μM PS/PO → 5 μM, 50 μL in 1 mL PBS

2 h incubation, wash, incubate 15 min Hoechst, wash, image

20x objective water flow, 6 images per well

Analysis: whole cell + endosomes detection

Continued to page:

SIGNATURE

DATE

18/11/19

WITNESSED AND UNDERSTOOD BY

DATE