

# Generation of Knockout-Overexpression cell lines

ID: SP0002-U

<b>Authors</b>	<sup>1</sup> Christoph Klimek, <sup>1</sup> Svenja Onstein, <sup>1</sup> Barbara Barbosa, <sup>1</sup> Alvaro Ingles-Prieto, <sup>1</sup> Anna Skucha, <sup>1</sup> Gernot Wolf
<b>Affiliations</b>	<sup>1</sup> CeMM Research Center for Molecular Medicine. Vienna, AT

## Protocol description

Protocol for the generation of knockout (KO) overexpression (OE) cell line pools via lentiviral transduction.

## Materials

<b>Biological materials:</b>	<ul style="list-style-type: none"><li>▪ Envelope plasmid pMD2.G (Addgene, 12259)</li><li>▪ Packaging plasmid psPax2 (Addgene, 12260)</li><li>▪ Lentiviral plasmid expressing HA-Strep-tagged codon-optimized cDNA (modified from pCW57.1, Addgene #41393)</li><li>▪ HEK293T cells</li><li>▪ Knockout cell line to be transduced</li></ul>
<b>Reagents:</b>	<ul style="list-style-type: none"><li>▪ Lipofectamine 3000 (ThermoFisher, L3000001)</li><li>▪ Blasticidin (Invivogen, ant-bl-5b)</li><li>▪ Opti-MEM (ThermoFisher, 31985070)</li><li>▪ Polybrene (Sigma-Aldrich, H9268)</li></ul>
<b>Equipment:</b>	<ul style="list-style-type: none"><li>▪ Sterile tissue culture hood, general lab equipment</li></ul>

## Reagents setup

### 5 mg/ml stock solution of Polybrene

Dissolve 50 mg of Polybrene in 10 ml of ddH<sub>2</sub>O and sterilize through a 0.22- $\mu$ m filter. Make 500  $\mu$ l aliquots of the solution and store at -20°C.

## Procedure

### ▪ Transfection - Day 1

1. Seed HEK293T cells in 100  $\mu$ l DMEM supplemented with FBS at a density of 15000 cells/well on a 96 well plate

### ▪ Transfection - Day 2

2. Check HEK293T cells (should be 50-70% confluent)
3. Prepare transfection reagents:
4. Lipofectamine MM (scale up according to planned transfections, add 10% extra)

	$\mu$ l per reaction
Optimem	7.15
Lipofectamine 3000	0.35

Total volume	7.5
--------------	-----

5. Let the Lipofectamine MM incubate for 10 min at RT
6. DNA MasterMix (MM) (scale up according to planned transfections, add 10% extra)

	µl per reaction
Optimem	4.8
pMD2.G (100 ng/ul)	0.3
psPax2 (100 ng/ul)	0.5
P3000 Reagent	0.3
Total volume	6.0

7. Pipette 6 µl of DNA MM per reaction into clean PCR strips or plates
8. Add 1.5 µl of overexpression plasmid (concentration: 50-100 ng/µl) to 6 µl of DNA MM
9. Mix 6 µl of Lipofectamine MM per reaction with each DNA MM
10. Incubate at RT for 5-15 minutes
11. Add 10 µl of mix to HEK293T cells (keep pipette tips in medium while dispensing but do not touch the bottom of the wells)

▪ **Transfection - Day 3**

12. Add 120 µl fresh medium to transfected HEK293T cells (pipette very slowly to avoid detachment of HEK293T cells)

▪ **Transfection - Day 4**

13. Transfer 200 µl of the virus-containing medium to a safe-lock tube and immediately freeze at -80°C

▪ **Transduction - Day 1**

1. Seed the cell line to be transduced at  $1.5 \times 10^4$  cells/well in a 96 well plate (Note: add one more well than needed as a negative control)

▪ **Transduction - Day 2**

2. Thaw the safe-lock tube with the virus-containing medium at RT

3. Change the medium of the cell line to 100 µl of medium with 10 µg/ml polybrene
  4. Add 100 µl of thawed lentivirus to the cells
- **Transduction - Day 3**
    5. If transduced cells are confluent: split transduced cells and transfer 200 µl to a 24 well plate
    6. If not confluent: change medium on the cells
  - **Transduction - Day 4**
    7. Change medium to selection medium (10 µg/ml Blasticidin)
    8. If cells are still on a 96 well plate, transfer to a 24 well plate and start Blasticidin selection the next day
  - **Transduction - Day 5-12**
    9. Monitor cells and split when reaching 80-100%. Keep under Blasticidin selection. Change medium every 3-4 days depending on confluency
  - **Transduction - Day 13-15**
    1. Observe negative control (untransduced cells), remove Blasticidin from transduced cells if all cells in neg. control are dead
    2. Monitor cells and expand when reaching 80-100% confluency (after at least 10 days of selection)
    3. Expand cells to 10 or 15 cm dishes (depending on how many freeze stocks are needed).

## **Additional notes**

## **Data/reagent availability**

<https://re-solute.eu/resources/reagents>

## **References**

Please write to [contact@re-solute.eu](mailto:contact@re-solute.eu) in case of questions or errors.