

Near *in vivo* microfluidic system to analyze Chemicals and NanoMaterials (CNMs) metabolic response

Prof. Maria Cristina Albertini, PhD

Associate Professor of General and Clinical Pathology


Department of Biomolecular Sciences (DISB)

University of Urbino Carlo Bo, Urbino (PU), Italy.

E-mail: maria.albertini@uniurb.it



New Approach Methodologies for risk assessment



- _CNMs
- _Scientific disciplines
- _Studies and analyses performed
- _Computer science and technology exploitation of primary research results



Lodz University of Technology



1506
UNIVERSITÀ
DEGLI STUDI
DI URBINO
CARLO BO



Funded by
the European Union



CheMatSustain has the support of the representatives of JRC of the European Commission/EURL ECVAM (European Reference Laboratory for Alternatives to animal testing).
<https://www.youtube.com/watch?v=LFayrdVGFIY>

- EU Reference Laboratory for alternatives to animal testing
- Promotes and facilitates the use of non-animal methods in research




Microfluidic devices can be used to create **models of blood circulation**:

- mimicking the behavior (shear stress) of the human circulatory system in a lab setting
- to study blood flow and cell behavior in a controlled environment
- to study cytotoxicity and inflammation processes


Static



Microfluidic



Blood flow



Shear stress: is the force per unit area that is created when a tangential force (blood flow) acts on a surface (endothelium).


Endothelium is highly sensitive to hemodynamic shear stresses that act at the vessel luminal surface in the direction of blood flow. **Physiological shear stress has been demonstrated to maintain endothelial quiescence and integrity.**

Experimental conditions:

9-13 June 2025

Depending on the aim of the project (toxicity analyses?) different cell types; time treatments; metabolic response have to be considered to set-up the microfluic system

Example: CheMatSustain



EA.hy926 cells: established by fusing primary human umbilical vein cells with a thioguanine-resistant clone of A549

Seeded cells	Time	Cell number	Cell death (Trypan blue)
60×10^3	3 days	$\pm 500 \times 10^3$	$\pm 4\%$

Our experiments last over 4 consecutive days (plus one washing day dedicated to all the components of the system).


MultiDyn (MD)



Chamber ($1,8 \text{ cm}^2$)

Fibronectin ($1 \mu\text{g}/\text{cm}^2$)


10x Magnification




Strengthening Research: UNIURB and TUL Collaborate on CNM Toxicity Testing

Toxicity assessment of nanoparticles through *in vitro* and *near in vivo* tests, employing human endothelial cells (EA.hy926 cells):

- Cytotoxicity; transcriptomics; proteomics
- inflammatory response (TLR-4 signaling by miR-146a; IRAK-1 and IL-6)



Lodz University of Technology



Day 1

Coating MD chambers

Seeding 6.0×10^4 EA.hy926 cells in 0,5 ml DMEM medium

Day 2

Flow activation at the rate of $250 \mu\text{l}/\text{min}$ for 24 h (conditioning)

Day 3

CNMs treatment for 24 h

Day 4

Cell detachment, cell viability and storing (-80°C)

Day 5

- Washing and sterilization of the microfluidic system (70% ethanol + UV)



Thank you!



<https://chematsustain.eu/>



<https://www.ivtech.it/>

The project CheMatSustain under No. 101137990 has received funding from the European Union under the Horizon Europe Programme.

