

Mass Spectrum SmartFormula Report

Analysis Info

Analysis Name E:\Data2\Kilian\KC-5-bpinPhCy0000001.d Acquisition Date 2017-03-23 18:18:59
Method Tune_low_pos.m Operator Carin Larsson
Sample Name PhCyBpin Instrument / Ser# microTOF 125
Comment kc-5-062

Acquisition Parameter

Source Type ESI
Focus Not active
Scan Begin 50 m/z
Scan End 1000 m/z
Ion Polarity Positive
Set Nebulizer 0.3 Bar
Set Dry Heater 180 °C
Set Dry Gas 4.0 l/min
Set Divert Valve Source

