

Federico Pazzaglia, Dr. Wouter Van Lysebettens, Xavier De Becker, Prof. Jeroen Missinne, Prof. Robrecht Raedt, Prof. Sandra Van Vlierberghe

CENTRE OF MICROSYSTEMS TECHNOLOGIES (CMST)



INTRODUCTION

Optogenetics (OG) and optopharmacology (OP) are powerful preclinical approaches relying on the illumination of tissues:


- Study of neuron networks behavior (OG)
- Activation and inhibition of genetically modified light responsive cells (eg. Neural networks) (OG)
- Selective drug release (OP)

OG and OP principles can be employed in closed-loop therapies for epilepsy where:


1. Electrodes record the EEG
2. A microcontroller processes the EEG
3. The optical probes will activate/inhibit the involved neurons (OG) or uncage drugs (OP)



INHIBITION OF EPILEPTIFORM AFTERDISCHARGES IN RAT BY BLUE LIGHT ILLUMINATION. BLUE LINE REPRESENTS THE ILLUMINATION ON [1]



NANOIMPRINTING LITHOGRAPHY PROCESS FLOW



HGs BENDING STIFFNESS VARIATION [2]

IMPLANTABILITY & BIOCOMPATIBILITY

Hydrogel (HG) biomaterials may be used as multifunctional coating:

- Soft and biocompatible coating
- Mechanical support during implantation

High stiffness is required to allow the probe penetration through tissues but foreign body reaction is limited by low stiffness probes.

Some HG (e.g PAAm-Alg, Acrylamide-Alginate) have swollen stiffness (kPa range) much lower than dehydrated state stiffness (Mpa range).

- Transition Time ~10'
- Demonstrated in agarose phantom brain

CHALLENGES

Chronic implantations of optical probes need to meet a number of requirements to guarantee a successful operation of the device.

1. Foreign body reaction (FBR) minimization
 - Shank Dimensions (<100µm)
 - Low stiffness
2. Optical Requirements
 - Low optical losses (<0.5dB/cm)
 - Irradiance up to 200mW/mm²
3. Surgical protocols compatibility
4. Limited heating of tissues during stimulation
5. Light source Integration


FABRICATION

Nano Imprinting Lithography (NIL)


- Pattern replication via Molding
- Inexpensive
- Fast prototyping

Process overview:


- Master Mold fabrication
 - Better control during imprinting
 - Lower probability to damage the master mold
- Imprint the pattern on Soft Stamp
 - Deposition of the optical material
 - Patterning via soft stamp imprinting
 - UV polymerization of the optical material



FABRICATED WAVEGUIDES TRANSMITTING LIGHT



SEM MICROGRAPH OF SOFT STAMPS



CROSS SECTION OF THE NIL FABRICATED WAVEGUIDES

CONCLUSIONS

- NIL technique is capable to produce flexible waveguides with overall dimension (thickness and width) smaller than 100µm
- Optical measurements reported more than 60% of light is transmitted (for WGs wider than 5µm)
 - An efficient coupling will lead to lower losses
- Microscope images show low surface roughness of the side walls
 - Surface roughness is one main responsible for optical losses
- Biomaterial coating is currently under development



WANT TO KNOW MORE?
SCAN THE QR CODE!

Contact

Federico.pazzaglia@ugent.be



FACULTY OF ENGINEERING
AND ARCHITECTURE

References

- [1] Acharya, AR et al. *CNS Neurosci Ther.* 2023; 29: 907-916. doi: [10.1111/cns.14029](https://doi.org/10.1111/cns.14029)
[2] Park, S. et al. *Nat Commun* 12, 3435 (2021). <https://doi.org/10.1038/s41467-021-23802-9>

SCAN ME