A miniaturized and low-cost sub-nanosecond fluorescence lifetime detector based on and array of CMOS SPAD detectors

Joan Canals, Nil Franch, Joel Diéguez, V. Moro, O. Alonso, Anna Vilà, Angel Diéguez

Department of Electronic and Biomedical Engineering. University of Barcelona. Spain adieguez@el.ub.edu



Outline

1. Motivation

- 2. Fluorescence as a technique for diagnosis
- 3. Lifetime detection based on CMOS SPADs
- 4. Miniaturized and low-cost sub-ns device
- 5. Conclusions



• Healthcare is under pressure in industrialized countries





3 Fast diagnosis at the point of care

Cardiovascular diseases (like acute myocardial infarction): leading cause of death

Strokes: Fast diagnosis between ischemic (TPA effective if less than 4 h) or hemorrhagic

Sepsis: One in three patients who die in a hospital have sepsis





World Congress on Biosensors 2018, 12-15 June 2018 | Miami, Florida, USA

Motivation

Fluorescence is a optical transduction technique

Fluorescence

2 Is very sensitive and selective (because few materials emit and thanks to labelling)

3 It can be measured Intensity, Lifetime and FRET (distance between fluorescent dyes)





1 Time Resolved Fluorescence is characteristic of each fluorescent molecule and the chemical

composition of the environment (can be used to follow reaction mechanisms)

2 It does not depend on excitation light intensity and concentration of fluorophores

3 Commonly measured by **TCSPC** (Time Correlated Single Photon Counting)





World Congress on **Biosensors 2018**, 12-15 June 2018 | Miami, Florida, USA

Fluorescence

O CMOS SPADs work in Geiger mode (bias 10-20% above breakdown) and have 'infinite' gain

2 Response not proportional to amount of light

CMOS SPADs

3 Dark Count rate is an issue (~7000 cps in this work)

Organization Probalility in visible (~350-700 nm)





UNIVERSITAT DE BARCELONA World Congress on **Biosensors 2018,** 12-15 June 2018 | Miami, Florida, USA



UNIVERSITAT DE BARCELONA



World Congress on **Biosensors 2018**, 12-15 June 2018 | Miami, Florida, USA

BARCELONA













BARCELONA



UNIVERSITAT DE BARCELONA



5000 excitation pulses t_{on}/t_{off} 30/20 ns total measurement time of 250 μ s

Minimum concentration $1/8 \mu M$ (low signal-to-noise ratio impeded to measure below)

Lifetime of τ =36.1±0.105 ns Agrees quite well with the reported value of 34.6 ns.

$$\tau = \frac{N(\Sigma t_i^2) - (\Sigma t_i)^2}{N \Sigma t_i \ln I_i - (\Sigma t_i)(\ln I_i)}$$

Q. Chen, A. Kiraz, X. Fan, Optofluidic FRET lasers using aqueous quantum dots as donors, Supplementary material, Lab. Chip., vol. 12, no. 2, pp. 353-359, 2016.



Conclusions

An array of 5 CMOS SPADs has been used to build a miniaturized fluorescence lifetime detector that could be used as a POC device for time resolved fluorescence

Is able to detect with >20% efficiency photons in the visible

2 The sample is loaded through a set of microfluidic channels
This could enable the system to perform multiplexed measurements with several biomarkers/labels
3 The system was tested with Quantum Dots in 20nL, down to concentration of 125nM.

Better resolution could be obtained implementing the SPADs in a process with lower noise





We acknowledge the Catalan government by financing this research under the contract with ACC1Ó 'Portable device for molecular diagnosis' coded **VALTEC13-1-0020-00**, and the contract with AGAUR 'Chip for molecular diagnosis' coded **2014LLAVOR00003**. We also acknowledge the Spanish government by financing this research with project coded **FPA2015-71292-C2-2-P**. We acknowledge the European Union by contract **737089**.





