MicroNanoFluidics 2018

15-16 March, Grenoble

On-chip Flow Cell Sorting System Based on High frequency Dielectrophoresis implemented on CMOS technology

Rémi MANCZAK - PostDoc (XLIM) Arnaud POTHIER (XLIM) Claire DALMAY, Pierre BLONDY (XLIM) Fabrice LALLOUE, Sofiane SAADA - PostDoc (HCP)








CSC location in tumor

Normal brain tissue Normal blood vessels












➡ How create CSC grow conditions in vitro ?

Adjusting culture media and removing 0_{2:} Stringent grow conditions










Phenotypical characterization of cell-pools





Sumcastec 🔰






Sensing module objectives

Recognize (distinguish from others) targeted cells Characterizing physical properties (<u>dielectrical</u>, mechanical, densitity, deformability...) of dissociated cells set in suspension

Isolate targeted cells (at least remove a maximum of unwanted cells)











Approach investigated

Combined hydro-fluidic & ElectroMagnetic (EM) manipulation



Cell are dynamically sorted depending their "susceptibility" to <u>specific</u> EM signal



HCP

Homéostasie Cellulaire et Pathologies Dielectrophoresis **(DEP)** approach








DEP - Methodology and measurement system






Sumcastec Xlin

- pathological state
- aggressiveness level


HCP Homéostasie Cellulaire et Pathologies

DEP - Methodology and measurement system



Frequency range variation – **UHF frequencies**

8












Crossover frequencies study

<u>Cell line</u> : Glioblastoma Ln18 (derived from malignant gliomas adult patients)







Crossover frequencies study


• **Standard** medium 95 cells analysed

• **Defined** medium 81 cells analysed

Note:














10

11

UHF \rightarrow Intracellular properties dominate the cell response








HCP Homéostasie Cellulaire et Pathologies



Combined hydro-fluidic & ElectroMagnetic DEP cell manipulation








HCP

Homéostasie Cellulaire et Pathologies




Combined hydro-fluidic & ElectroMagnetic DEP cell manipulation

First prototypes based on IHP technology with PDMS cap for microfluidic



\sim 9 μm thick electrodes embedded in the $\mu fluidic$ channel

13




Fabricated in IHP clean room




First trials with flowing cells

Effect of applied DEP E field on particle trajectory

High efficiency on low cell velocity motion (< 0.1mm/s)





Stacked image of cell trajectory deviation

Particle is quickly repealed far the high intensity E field region to finally travel in low intensity one










Optimization of cell injection inside the cytometer

15






Microfluidic








Preliminary Experiments on LN18 cells (standard medium)



F_{DEP}< Fcross-over As expected cells are distributed on the edge of the channel attracted and trapped by high intensity E field areas












Preliminary Experiments on LN18 cells (standard medium)

Sumcastec 🗙










Particle sorting capability



Evaluation of sensor sorting capability to separate different type of particles

- With a 50MHz DEP signal: all polystyrene beads react in Neg DEP
 - most of LN18 cells react in Pos DEP
 - very few cells seem not be deviated (dead or damaged cells?)



Stacked image of LN18 cells and polystyrene beads with a 50MHz DEP signal

At the sensor output particle distribution testifies of its good sorting ability








Conclusion

Sorting principle validated



> Cell injection condition and suitable cell speed still under optimization

- Future tests will imply new high speed particle tracking monitoring system
- Future Implementation of fluorescence imaging capability to follow labelled cells with die

> Sorting of cell mixtures according to their crossover frequencies







19

Acknowledgments

We would like to thank ours collaborators and especially SUMCASTEC project partners (FET OPEN Research and innovation project) grant agreement N° 737164.





20

Biology







