



PRIFYSGOL
BANGOR
UNIVERSITY



Electro-permeabilization of isolated cancer stem cells with a push-pull configuration of high power MOSFETs

Ilan Wyn Davies

Creo Medical / Bangor University

Introduction to the Project





Introduction to the Project



PRIFYSGOL
BANGOR
UNIVERSITY



Agenzia nazionale per le nuove tecnologie,
l'energia e lo sviluppo economico sostenibile



UNIVERSITÀ
DEGLI STUDI
DI PADOVA

- Part of European Union's Horizon 2020 research and innovation program: **Semiconductor based Ultrawideband Micromanipulation of Cancer STEM Cells or SUMCASTEC**
- SUMCASTEC explores a new approach for real time isolation and neutralization of Cancer Stem Cells (CSCs).
<http://www.sumcastec.eu>
- CSCs are associated with Glioblastoma Multiforme and Medulloblastoma.
- A project deliverable: to develop an off-chip pulsed Electric-Field (EF) generator for the Electroporation (EP) of CSCs.





Project Deliverables

- To deliver a generator capable of Electroporation
- Generate pulse of amplitude in excess of 1 kV, with pulse widths of 100 ns to 300ns to a 50Ω load.
- Minimisation of overshoot and ringing (Flat-Top pulses)
- Investigate various pulse parameters of SUMCASTECs Pulse Generator (SPG) on CSCs. (parameters as: pulse duration, repetition frequency, number of pulses)
- CSCs placed into a cuvette with a 50ohm buffer solution
- Develop an Non-thermal EP approach


Background






What are Cancer Stem Cells?

- Tumour growth has been linked to the presence of CSCs
- Main embodiment of a tumour with cancer cells
- Resistance to Chemotherapy (general chemotherapeutic agents) & Radiation Therapy
- Posses the ability of recovering from cytotoxic therapy and slowly repopulate the tumour, leading to relapse.




Tumour = Weed





What is Electroporation?



- Alternative physical technique for non-thermal ablation
- Use precisely controlled high amplitude **pulsed electric fields** of **short duration** to alter a cell's transmembrane potential.
- Results in permeabilizing the cell's plasma membrane and disturbing intercellular homeostasis.
- The resultant permeabilization of cell plasma membrane can be reversible or irreversible

SUMCASTEC EP Pulse Generator (SPG)






Instrumentation

- Push-pull switching of High Voltage, fast switching MOSFETs
- Driven by opto-isolators with comparable switching times
- High enough current to charge up gate-source and gate drain capacitances

$$i = C \frac{dV}{dt}$$


- High-side (HS) MOSFET determine pulse width. Low-side (LS) MOSFET complimentary of High-side MOSFET





Characterization

- Flat pulses free from ringing and overshoot
- Developed EP generator performance exceeds the LTSpice simulation



- Increase of gate voltage from 15 V to 20 V results in increased pulse amplitude.
- Pulse amplitude is unaffected throughout its operating repetition frequencies (1-50 Hz)
- Optimized for pulses of 100 ns to 300 ns in width, and amplitudes in excess of 1kV

Project Setup



Instruments

SUMCASTEC Pulse Generator (SPG)



Pulse repetition Frequency: 1 - 50 Hz
Pulse Width Generated: 80 ns - 1 μ s*


Number of pulses generated: 1 - ∞
infinite (continuous wave)

Pulse amplitudes: 280V - 1100V**

N-connector and banana sockets output.

*(increment of 10ns between 80 - 400 ns and 20ns increment 400ns - 1 μ s) ** With 50 Ω Load Designed and Build by Creo Medical

ENEA Cuvette Housing Unit*




Exposes CSCs in cuvette to EP
Compatible with CSCs and standard N-connectors

Allow real time monitoring of EP HV pulse

*Designed and constructed by A. Zambotti and ENEA

ENEA 50 Ω artificial buffer solution



Butten (4)PBS + H₂O + Sucrose [1]
Proportions for 100 mL of 0.3 S/m (50 Ω) buffer solution:


- 20 mL of PBS (phosphate saline buffer)
- 80 mL of H₂O (distilled water)
- 8.2 g of sucrose (to balance out Osmotic Pressure with the cells)

Made by ENEA



Dynamic monitoring of pulse
during cell EP or IRE

Bench Study Set-up



**SUMCASTEC
Pulse
Generator
(SPG)**

**ENEA Cuvette
Housing Unit**


*Bio-rads 0.1cm gap Cuvette
containing CSCs and ENEA
50Ω artificial buffer solution*

Results and Analysis






Pulsed Electric-Field on CSC



- AI ENEA artificial 50Ω is comparable to waveform measured with a 50Ω resistor.
- 100 ns, 200ns and 300 ns pulse waveforms measured across the EP cuvette containing CSCs suspended in 50Ω , 0.3 S/m buffer solution



Pulsed Electric-Field on CSC




- 300ns pulses with various buffer solution at load
- Shape of the pulse is non-affected
- Demonstrating broadband matching performance




Permeabilization Test on CSCs

Permeabilized





Unpermeabilized



- 1. CSCs placed in buffer solution with YOPRO-1, a fluorescent dye (green in colour)
- 2. CSCs exposed to EP
- 3. Exam CSCs under fluorescent light:
 - a. Permeabilized CSCs green
 - b. Unpermeabilized: CSCs transparent



Permeabilization Test Result



300 ns, 1.2 MV/m, 20 Pulses
1 Hz = 90 % permeabilized



100 ns, 1.0 MV/m, 20 Pulses
0% permeabilized



Non-thermal Effect

Pulse Width, P_w (ns)	Amplitude V (kV)	Repetition Frequency, f (Hz)	Energy, E (mJ)	Temperature Change, ΔT ($\mu^\circ C$)
100	1.0	1	2.00	4.8
100	1.0	50	100.00	239.2
200	1.2	1	5.76	13.8
200	1.2	50	288.00	689.0
300	1.2	1	8.64	20.7
300	1.2	50	432.00	1033.0

Non-thermal effect. of $1.0 \times 10^{-3} \text{ }^\circ C$ ($100 \mu^\circ C$)

- D is duty cycle (ratio)
- E is energy (J)
- C is heat coefficient, $4.18 \text{ J/g/}^\circ C$, same as water
- L (indicating volume) is millilitres (the cuvette can hold 0.1 mL of solution)

$$\Delta T = \frac{\left(\frac{V^2}{Z}\right) \cdot P_w \cdot D}{C \cdot L}$$

Discussion and Conclusion





Project Deliverables

- ✓ Successful permeabilization of the CSCs
 - ✓ Real-time pulse visualization
 - ✓ Positive results obtained in matching strategy
 - ✓ Non-thermal permeabilization of CSCs
-
- Results do not reflect whether the CSCs are dead or alive
 - Unknown if EP was reversible or irreversible
 - CSCs exposed to 100 ns pulses could have been permeabilized
 - Viability studies are required to complement these preliminary experiments in the future



Acknowledgment:

The project received funding from the European Union's Horizon 2020 research and innovation program under grant agreement No.737164. The authors would like to thank partners of the SUMCASTEC project for the successful collaboration. We thank A. Casciati, M. Tanori, B. Tanno and M. Macuso from ENEA for the support in cell culturing and microscopy.

Referneces:

[1] C. Merla et al. "SUMCASTEC_180123_NA_protocolWP3_protocol_.pdf_Rome_C.M.Merla_Partners and public_NA", Zenodo, 2018.

For more information visit the SUMCASTEC website:

www.sumcastec.eu



Our project SUMCASTEC
was made possible
thanks to #H2020 funding

ENEA

Agenzia nazionale per le nuove tecnologie,
l'energia e lo sviluppo economico sostenibile

€30 billion is still
available in the 2018-20
Work Programme!

#InvestEUresearch



www.sumcastec.eu



UNIVERSITÀ
DEGLI STUDI
DI PADOVA