

# BioRECO<sub>2</sub>VER

8th Conference on CO<sub>2</sub> as Feedstock for Fuels, Chemistry and Polymers

Heleen De Wever and project partners, 24 March 2020




Horizon 2020  
European Union Funding  
for Research & Innovation

This project has received funding from the European Union's Horizon 2020 research and innovation programme under grant agreement No. 760431.

## General information

- Title Biological routes for CO<sub>2</sub> conversion into chemical building blocks
- Acronym BioRECO<sub>2</sub>VER
- Work program topic BIOTEC-5-2017: Microbial platforms for CO<sub>2</sub>-reuse processes
- Type of action Research and Innovation Action
- Start date 1<sup>st</sup> January 2018
- End date 31<sup>th</sup> December 2021
- EU budget € 6,812,187.50
- Coordinator VITO

# Overall project concept




# Project Consortium

- 4 Industries
- 4 SMEs
- 2 universities
- 2 RTOs



# CO<sub>2</sub> capture and pretreatment - An overview



# Enzyme improvement by directed evolution

## Library construction:

1000 epCA variants generated with error-prone PCR


## Library screening:

- Primary screening to select active variants
- Secondary screening to select variants with resistance to inhibitors

Scaled-up production of most promising variants


Sequencing for identification of mutations

- **3 mutants showed 50% increased resistance to flue gas inhibitors**



# Enzyme improvement by immobilization

Applied techniques for CA immobilization:








# Development of a solvent with competitive absorption and desorption properties

Screening of different AAIL:tertiary amine blends in small scale resulted in selection of solvent with:

- **5-fold higher initial rate**
- **2-fold higher CO<sub>2</sub> load compared to MDEA**
- **2-fold higher regeneration at 80°C**
- **>15% reduction in desorption T compared to MEA**



MDEA: Methyl diethanolamine; MEA: Monoethanolamine



## CO<sub>2</sub> capture in large-scale packed bed absorption equipment

Scaled-up trials revealed even higher benefit with use of the developed solvent:

Solvent	Relative $K_G a$ (%)
25% MEA	100
25% MDEA	3.07
25% AAIL: amine blend	30.86



- **10-fold increase in mass transfer coefficient ( $K_G A$ ) compared to MDEA**



1m (80mm ID) packed column (Raschig rings)

MDEA: Methyl diethanolamine; MEA: Monoethanolamine

# Integration of CA with novel hybrid solvent for efficient CO<sub>2</sub> capture



- Reduced operation times by 25%
- 32% increase in captured CO<sub>2</sub> compared to the non-enzymatic reaction

MDEA: Methyl diethanolamine

## CO<sub>2</sub> capture and pretreatment

### Conclusions:

- An ultrastable CA was improved by protein engineering and immobilization for increasing stability towards harsh and high temperature environment
- An enzyme compatible novel hybrid solvent was developed with competitive absorption and desorption properties

### Next steps:


- Integration of *immobilized improved epCA variants* with the developed *novel hybrid solvent* for *large scale absorption* using *flue gas* resulting in a concentrated gas stream of at *least 95% CO<sub>2</sub>*

# Microbial CO<sub>2</sub> conversion

- 3 microbial platforms

Microbial platforms		T range	O <sub>2</sub> tolerance	Target product	Partner
Autotrophic	Clostridial strain	Mesophilic	Anaerobic	Isobutene	GLOBAL BIOENERGIES
	<i>Cupriavidus necator</i>	Mesophilic	Aerobic	Lactate	
Capnophilic	<i>Thermotoga neapolitana</i>	Hyper-thermophilic	Strictly anaerobic	Lactate + H <sub>2</sub>	 Consiglio Nazionale delle Ricerche


# Wood-Ljungdahl and Isobutene Pathway



- WLP is the most efficient carbon fixation pathway
- Isobutene is derived from central metabolic precursor acetyl-CoA
- Isobutene pathway is Redox-neutral and does not require directly ATP

## Isobutene production under autotrophic conditions

- Successful implementation of Isobutene pathway proven on protein level
- Isobutene production from CO, CO<sub>2</sub> and H<sub>2</sub> mixture and from CO<sub>2</sub>/H<sub>2</sub> mixture
- Large improvement of production during the BioRECO<sub>2</sub>VER project



# Aerobic lactate production

Metabolic engineering of *Cupriavidus necator* strain to improve lactate production from CO<sub>2</sub>

## 1. Improve lactate dehydrogenase (LDH) activity


- LDH activity increased by a factor 180
- Specific lactate production rate increased by a factor 4

## 2. Increase pyruvate availability by deletion of competitive pathways

- Different competitive pathways deletions evaluated

## 3. Block lactate re-consumption


- Unexpected pathway identified by transcriptomic study
- Lactate re-consumption issue solved by deletion work
- Patent application filed



Lithoautotrophic metabolism



# Capnophilic lactate production




**Capnophilic Lactic Fermentation (CLF) pathway:**  
*Thermotoga neapolitana*-based platform to gain value from CO<sub>2</sub> and waste by production of L-lactate & H<sub>2</sub>

- Newly discovered pathway (ChemSusChem, 2014, 7, 2678-2683)
- Dissection of anabolic branch of CLF
- Increase metabolic flow from CO<sub>2</sub> & acetate to lactate
- Proof of concept net CO<sub>2</sub> fixation in lactic acid
- Explore feeding of exogenous acetate as C2 unit

## Capnophilic lactate production

- Selection of two model strains with increased CLF productivity
  - *T. neapolitana* subsp. *capnolactica* (DSMZ33003) is a mutant of DSMZ 4359<sup>T</sup> strain and has been generated in our laboratory under saturating concentrations of CO<sub>2</sub>
  - **RQ7** is a *Thermotoga* sp. strain isolated from hot sea-floor Ribeira Quente (the Azores)

➔ Complete genome sequences are available






\*100 % more than reference strain DSMZ 4359<sup>T</sup>

# Capnophilic lactate production

## DSM33003



- heterologous expression AcCoA synthetase (*acs*) of *Thermus thermophilus* (mutant **CL-11-AS**)

CLF performances under standard medium + 20 mM acetate



- Comparable growth curves
- Same glucose consumption rate
- Decrease H<sub>2</sub> production
- 100% increase lactate production


# Capnophilic lactate production




Increase of CO<sub>2</sub>-  
derived lactate of  
70%

## UPTAKE

### Wild type



### CL-11-AS




# Microbial CO<sub>2</sub> conversion


- 2 technologies




Bioelectrochemical systems



[High pressure] fermentors




# Bioelectrochemical systems (BES)












## Cell configurations

**Clostridium platform:**  
**Double chamber**  
 - Anaerobic conditions

**Cupriavidus platform:**  
**Single chamber**  
 - Membraneless system  
 - Coexistence of  $O_2$  and  $H_2$




Maximize *in situ* H<sub>2</sub> production






~~O<sub>2</sub>~~ 50°C

# Thermophilic process




Thermophilic MEC



High T°

- Higher reaction rates
- Less risk of contamination
- More product specificity



# Thermophilic process

Set-up of mild thermophilic systems

## Chronology

Reactors 1 and 2 (280 d operation)



After 70 d, inoculation of


Reactors 3 and 4 (210 d operation)



After 160 d, inoculation of

Reactors 5 and 6 (50 d operation)


- Max production rate: 28 g acetate m<sup>-2</sup> d<sup>-1</sup>
- Coulombic efficiency: 80-90%



## Aerobic process


- *Cupriavidus* lactate fermentation in a single chamber BES
  - 2-stage process
  - Expected stoichiometries

### Biomass production *Cupriavidus*




Gas mix CO<sub>2</sub>/O<sub>2</sub>/H<sub>2</sub> = 16/20/64 (literature)

### Lactate production




# Aerobic process

- Growth of wild type *Cupriavidus* strain
  - Stainless steel electrodes result in OD increase while graphite ones do not  
(SS: higher  $H_2$  evolution)  
(Graphite: sorption - peroxide)
  - Too negative voltage results in slower growth  
(foam formation)




## Aerobic process

- Tests with modified *Cupriavidus* strain
  - Conditions
    - Stainless steel electrodes
    - Avoid O<sub>2</sub> inhibition enzymes
    - Apply different voltages
    - Test different electrode surface/liquid ratio
  - Conclusions
    - Faster growth in reactors with higher electrode surface area
    - Fixed ratio O<sub>2</sub>/H<sub>2</sub> from water electrolysis has risk O<sub>2</sub> inhibition or H<sub>2</sub> limitations
    - Headspace composition to be optimized




# Bioelectrochemistry

- Next steps
  - Co-culture tests



- Benchmark with proprietary gas diffusion electrode



# High pressure fermentation

- Features
  - Operating pressure = 1 to 10 bara
  - H<sub>2</sub> and CO<sub>2</sub> as process gas; also for gas mixtures with O<sub>2</sub>, CH<sub>4</sub> and CO, real offgases
  - Online GC to monitor headspace composition
  - Separate electrical cabinet controlling battery of mass flow controllers for gas addition
  - Design allows 100% gas consumption
  - Integrated membrane filtration unit for cell retention
  - *In situ* product recovery
- ATEX compliance
- Fully automated – Labview
- Commissioned December 2019 – January 2020





Online GC




Fermentor skid



# High pressure fermentation

- Open-access review paper (see project website)
  - Different types of reactors discussed
  - H<sub>2</sub> and CO sensors in liquid: troublesome
  - Titrers and productivities need to be improved significantly to allow industrialization
  - Gas to ethanol most advanced
  - More optimization work required to produce non-native chemicals from C1 substrate gases



Contents lists available at [ScienceDirect](https://www.sciencedirect.com)

Bioresource Technology

journal homepage: [www.elsevier.com/locate/biortech](http://www.elsevier.com/locate/biortech)

Review

Effects of moderately elevated pressure on gas fermentation processes

Wouter Van Hecke<sup>a</sup>, Richard Bockrath<sup>b</sup>, Heleen De Wever<sup>a,\*</sup>

<sup>a</sup> Separation and Conversion Technology, Flemish Institute for Technological Research (VITO), Boeretang 200, Mol, Belgium

<sup>b</sup> Global Bioenergies, 5 rue Henri Desbruères, 91000 Evry, France


# Microbial conversions

## • Conclusions

- Established manipulation and genetic toolbox for *C. ljungdahlii* and implemented the isobutene pathway
- Isobutene production under autotrophic conditions
- Lactate re-consumption issue solved within *Cupriavidus* and lactate production improved by lactate dehydrogenase overexpression
- Selection of two model strains of *Thermotoga neapolitana*, DSM33003 and RQ7
  - productivity (increase lactic acid molarity)
  - genetic tools (transformable strains and amelioration of target steps)
- Test work with bioelectrochemical and high pressure set-ups initiated
- Further optimizations to improve titers and productivity

# Acknowledgements


*This project has received funding from the European Union's Horizon 2020 research and innovation programme under grant agreement No. 760431.*



# Thank you for your attention!

Contact:

heleen.deweever@vito.be



## Publications

- W. Van Hecke; R. Bockrath; H. De Wever (2019); *Effects of moderately elevated pressure on gas fermentation processes*, DOI: 10.1016/j.biortech.2019.122129
- V. Luongo; A. Palma; E. R. Rene; A. Fontana; F. Pirozzi; G. Esposito; P. N.L. Lens (2018); *Lactic acid recovery from a model of Thermotoga neapolitana fermentation broth using ion exchange resins in batch and fixed-bed reactors*, DOI:10.1080/01496395.2018.1520727
- G. Dreschke, G. d'Ippolito, A. Panico, P. N.L. Lens, G. Esposito, A. Fontana (2018); *Enhancement of hydrogen production rate by high biomass concentrations of Thermotoga neapolitana*, DOI: 10.5281/zenodo.3247830
- G. Nuzzo; S. Landi; E. Nunzia; E. Manzo; A. Fontana; G. d'Ippolito (2019); *Capnophilic Lactic Fermentation from Thermotoga neapolitana: A Resourceful Pathway to Obtain Almost Enantiopure L-lactic Acid*, DOI: 10.3390/fermentation5020034
- N. Pradhan; G. d'Ippolito; L. Dipasquale; G. Esposito; A. Panico; P.N.L. Lens; A. Fontana (2019); *Simultaneous synthesis of lactic acid and hydrogen from sugars via capnophilic lactic fermentation by Thermotoga neapolitana cf capnolactica*, DOI: 10.5281/zenodo.3247821



Horizon 2020  
European Union Funding  
for Research & Innovation

BioRECO<sub>2</sub>VER

*The sole responsibility for the content of this dissemination and communication activity lies with the authors. It does not necessarily reflect the opinion of the European Union (EU) Horizon 2020. The EU is not responsible for any use that may be made of the information contained therein.*