

Optimizing molecular assays to support early detection of *Xylella fastidiosa* in host plants

Chrysavgi I. Reppa¹, Paraskevas E. Glynos¹,
Charikleia D. Karafila¹, Maria C. Holeva¹

¹Benaki Phytopathological Institute, Department of
Phytopathology, Laboratory of Bacteriology, Kifissia, Greece

INTRODUCTION & OBJECTIVES

- Efficient sampling approaches combined with reliable, fast testing procedures are crucial factors in plant surveillance for early detection of *Xylella fastidiosa* (Xf).
- In the frame of the European project 'XF-ACTORS', work was performed on developing effective sampling methods to support reliable detection of Xf, particularly when plants are asymptomatic.

The **objectives** of this work were:

- to evaluate the efficiency of two methods for DNA extraction from plant material (1g), employing two commercial kits: one silica-column-based and one magnetic-bead-based, in combination with automated DNA-extraction platforms.
- to test the performance of both above-mentioned DNA extraction methods in pooled (10g) plant samples.
- to evaluate the sensitivity of an isothermal DNA amplification and detection system applied on crude plant material (homogenate) without DNA purification, using a portable fluorometer (AmpliFire, Douglas Scientific).

MATERIALS & METHODS

Detection method: qPCR

Plant matrix:

Olea europaea, *Prunus ceracifera*, *Ficus carica*, *Lavandula* sp., *Nerium* sp., *Pelargonium* sp., *Pistacia* sp., *Juniperus* sp. and *Tagetes* sp. (i.e. hosts of high economic importance & wide distribution in Greece, or major role in Xf epidemiology): 1g or 10g

Bacterial inoculum:

Inactivated Xf cells (CFBP 8477) were added in the plant matrix (spiked samples) at 10^6 , 10^5 , 10^4 , 10^3 , 10^2 cfu/ml

Sample preparation

DNA extraction by:

- DNeasy Mericon™ Food Standard Protocol [EPP0 PM7/24(4)] using Qiacube platform (Qiagen)
- Maxwell® RSC PureFood GMO and Authentication Kit AS1600 using Maxwell® RSC platform (Promega)

DNA extracts from Xf spiked plant matrix were assessed in regard to their:

- quality (A_{260}/A_{280})
- quantity (ng/μl)
- performance in qPCR: Harper *et al.* (2010; erratum 2013); Francis *et al.* (2006)

Detection method: isothermal DNA amplification

Plant matrix:

Polygala myrtifolia (indicator host for monitoring the presence of insect vectors of Xf) and *Myrtus microphylla* (species with strong inhibition in conventional PCR)

Bacterial inoculum:

Inactivated Xf cells (CFBP 8477) were added in the plant matrix (spiked samples) at 5 to 1250 cells/reaction

Sample preparation


- Crude plant homogenate with no DNA extraction
- Kit of AmplifyRP XRT, using the portable fluorometer (AmpliFire, Douglas Scientific)

Conditions for testing:

- 'standard' (ST), recommended by manufacturer
- ST+incubation of homogenate in AMP1 buffer for 15 min
- ST+incubation of homogenate in AMP1 buffer for 60 min
- ST+low speed centrifugation of homogenate (5130 rpm/5 min)

Positive control: Xf cells (inactivated; 625 cells/reaction, i.e. 10^7 cfu/ml)



Negative control: healthy *Myrtus* sp. tissue





RESULTS

For spiked plant tissue of 1g

- Concentration: DNeasy Mericon < Maxwell RSC
A260/280: DNeasy Mericon > Maxwell RSC
- Cqs for target **COX** \approx Cqs for target **18S**, for *Olea*, for both kits
Cqs for target **COX** > Cqs for target **18S**, for *Juniperus*, *Pelargonium*, *Prunus*, *Tagetes*, for both kits
- Cqs by qPCR **Harper et al.** < Cqs by qPCR **Francis et al.** for all hosts, and both kits
- qPCR **Harper et al.**:
Cqs after DNeasy **Mericon** > Cqs after **Maxwell RSC**, for all hosts spiked at 10^6 to 10^4 cfu/ml. The opposite was true for *Pelargonium* 10^3 cfu/ml, *Olea* and *Juniperus* 10^2 cfu/ml



- qPCR **Francis et al.**:
Cqs after DNeasy **Mericon** varied among host plants and inoculum level from higher to lower compared to Cqs after **Maxwell RSC**



RESULTS

For spiked plant tissue of 10g


- Concentration: DNeasy Mericon < Maxwell RSC (especially *Lavandula* sp.)
A260/280: DNeasy Mericon > Maxwell RSC only for *Ficus* and *Lavandula*
- Cqs for target **COX** > Cqs for target **18S** (all hosts, both kits)
Cqs by qPCR **Harper et al.** < Cqs by qPCR **Francis et al.** (all hosts, both kits)
- qPCR **Harper et al.**:
Cqs after DNeasy **Mericon** > Cqs after **Maxwell** RSC, for *Ficus* and *Pistacia*
Cqs after DNeasy **Mericon** varied among host plants and inoculum level from higher to lower compared to Cqs after **Maxwell** RSC
- qPCR **Francis et al.**:
Cqs after **DNeasy Mericon** > Cqs after **Maxwell** RSC, for *Ficus* and *Pistacia*
Cqs =0 after Maxwell for *Lavandula*, while detection possible after DNeasy Mericon

MAIN CONCLUSIONS & ACKNOWLEDGEMENTS


- No single DNA extraction protocol can be recommended for detection of Xf in all plant species. The choice of method depends also on the amount of plant tissue to be analyzed.
- Harper et al. (2010; erratum 2013) qPCR protocol performed well with DNA extracts from large amount of tissue, which is significant for pooling samples together.
- Targeting eukaryotic 18S rDNA as internal control in qPCR is more efficient than targeting COX.
- Isothermal amplification system (AmpliFire) provided promising results for onsite Xf detection.

This work has received funding from the European Union's Horizon 2020 Research and Innovation Programme under the Grant Agreement N. 727987 "Xylella fastidiosa Active Containment Through a multidisciplinary-Oriented Research Strategy XF-ACTORS". Maxwell® RSC platform (Promega) and AmpliFire portable fluorometer (Agdia) were kindly provided by the manufacturers for a short testing period and we would like to express here our sincere thanks to them.

Isothermal DNA amplification



The sensitivity of the AmpliFire system was estimated at approximately 5×10^3 cfu/ml.



Modifications did not improve sensitivity. The best results were obtained with the use of the immunological cassette post amplification.