



DRAWN		Blood Pucker		
Jim	6/8/2020			
CHECKED		TITLE		
QA		All units in mm		
MFG		Material: Stainless Steel 304 or 316		
APPROVED		Remove all burrs and Break all edges		
		SIZE	DWG NO	REV
		A4	Blood pucker_3DModel	
		SCALE	SHEET 1 OF 1	
		1 : 1		

TITLE

PART NO

Blood puck_3DModel

REVISION

DESIGNER

Shop

ENGINEER

NOTES

2020 Update Blood Puck files to have a 3D model reflecting current device in use. And order more for rapid Phase in in the vosshall lab

THE INFORMATION AND/OR MATERIAL IN THIS DOCUMENT IS THE PROPERTY OF AND RESTRICTED INFORMATION AND/OR MATERIAL OF THE AUTHOR. THIS INFORMATION MAY NOT BE USED, REPRODUCED, PUBLISHED OR DISCLOSED TO OTHERS WITHOUT WRITTEN AUTHORIZATION. IT IS TO BE USED ONLY FOR MANUFACTURING ITEMS SPECIFIED WITHIN THE DOCUMENT.